

LYMPHOID NEOPLASIA

Standardized flow cytometry for highly sensitive MRD measurements in B-cell acute lymphoblastic leukemia

Prisca Theunissen,^{1,*} Ester Mejstrikova,^{2,*} Lukasz Sedek,³ Alita J. van der Sluijs-Gelling,⁴ Giuseppe Gaipa,⁵ Marius Bartels,⁶ Elaine Sobral da Costa,⁷ Michaela Kotrová,² Michaela Novakova,^{2,6} Edwin Sonneveld,⁴ Chiara Buracchi,⁵ Paola Bonaccorso,⁵ Elen Oliveira,⁷ Jeroen G. te Marvelde,¹ Tomasz Szczepanski,³ Ludovic Lhermitte,⁸ Ondrej Hrusak,² Quentin Leclercq,⁹ Georgiana Emilia Grigore,¹⁰ Eva Froňková,² Jan Trka,² Monika Brüggemann,⁶ Alberto Orfao,⁹ Jacques J. M. van Dongen,¹ and Vincent H. J. van der Velden,¹ on behalf of the EuroFlow Consortium

¹Department of Immunology, Erasmus MC, University Medical Center Rotterdam, Rotterdam, The Netherlands; ²Childhood Leukaemia Investigation Prague, Department of Paediatric Haematology and Oncology, Second Faculty of Medicine, Charles University, University Hospital Motol, Prague, Czech Republic; ³Department of Pediatric Hematology and Oncology, Zabrze, Medical University of Silesia, Katowice, Poland; ⁴Dutch Childhood Oncology Group, The Hague, The Netherlands; ⁵Centro Ricerca Tettamanti, Clinica Pediatrica Università di Milano Bicocca, Monza, Italy; ⁶Department of Hematology, University of Schleswig-Holstein, Campus Kiel, Kiel, Germany; ⁷Department of Pediatrics, Federal University of Rio de Janeiro, Rio de Janeiro, Brazil; ⁸Department of Hematology, Hôpital Necker-Enfants-Malades and Unité Mixte de Recherche Centre National de la Recherche Scientifique 8147, University of Paris Descartes, Paris, France; ⁹Cancer Research Center (Instituto de Biología Molecular y Celular del Cáncer–Consejo Superior de Investigaciones Científicas), Department of Medicine and Cytometry Service, University of Salamanca, Institute of Biomedical Research of Salamanca, Salamanca, Spain; and ¹⁰Cytognos SL, Salamanca, Spain

Key Points

- Standardized flow cytometry allows highly sensitive MRD measurements in virtually all BCP-ALL patients.
- If sufficient cells are measured (>4 million), flow cytometric MRD analysis is at least as sensitive as current PCR-based MRD methods.

A fully-standardized EuroFlow 8-color antibody panel and laboratory procedure was stepwise designed to measure minimal residual disease (MRD) in B-cell precursor (BCP) acute lymphoblastic leukemia (ALL) patients with a sensitivity of $\leq 10^{-5}$, comparable to real-time quantitative polymerase chain reaction (RQ-PCR)-based MRD detection via antigen-receptor rearrangements. Leukocyte markers and the corresponding antibodies and fluorochromes were selected based on their contribution in separating BCP-ALL cells from normal/regenerating BCP cells in multidimensional principal component analyses. After 5 multicenter design-test-evaluate-redesign phases with a total of 319 BCP-ALL patients at diagnosis, two 8-color antibody tubes were selected, which allowed separation between normal and malignant BCP cells in 99% of studied patients. These 2 tubes were tested with a new erythrocyte bulk-lysis protocol allowing acquisition of high cell numbers in 377 bone marrow follow-up samples of 178 BCP-ALL patients. Comparison with RQ-PCR-based MRD data showed a clear positive relation between the

percentage concordant cases and the number of cells acquired. For those samples with >4 million cells acquired, concordant results were obtained in 93% of samples. Most discordances were clarified upon high-throughput sequencing of antigen-receptor rearrangements and blind multicenter reanalysis of flow cytometric data, resulting in an unprecedented concordance of 98% (97% for samples with MRD < 0.01%). In conclusion, the fully standardized EuroFlow BCP-ALL MRD strategy is applicable in >98% of patients with sensitivities at least similar to RQ-PCR ($\leq 10^{-5}$), if sufficient cells ($>4 \times 10^6$, preferably more) are evaluated. (Blood. 2017;129(3):347-357)

Introduction

Most current treatment protocols for B-cell precursor (BCP) acute lymphoblastic leukemia (ALL) include minimal residual disease (MRD) measurements, generally based on polymerase chain reaction (PCR) analysis of rearranged antigen receptor genes.¹⁻³ Although flow cytometry (FCM) can be used for MRD detection as well,⁴⁻⁹ studies so far indicate that the specificity and sensitivity of FCM-MRD diagnostics are inferior to PCR-based MRD diagnostics.¹⁰⁻¹³ Nevertheless, we and others have recently shown that the use of 6- or 7-color immunostainings combined with the introduction of new markers and new marker combinations significantly improved FCM-MRD analysis

in BCP-ALL patients.^{10,12} These improvements were particularly related to specificity, whereas the sensitivity still appeared to be lower than for the PCR-based methods. To further improve FCM-based MRD diagnostics, more objective and efficient discrimination of BCP-ALL cells from normal BCP cells and improved sample preparation procedures for acquisition of larger numbers of cells are a prerequisite.

Eight-color immunostainings may contribute to improve flow cytometric MRD detection in BCP-ALL patients. Recently, an 8-color antibody tube was developed in the ALL-REZ-BFM 2002 trial.¹⁴ This tube contained 7 antibodies (CD10, CD19, CD20, CD22, CD34,

Submitted 5 July 2016; accepted 23 November 2016. Prepublished online as *Blood* First Edition paper, 30 November 2016; DOI 10.1182/blood-2016-07-726307.

*P.T. and E.M. contributed equally to this study.

The online version of this article contains a data supplement.

The publication costs of this article were defrayed in part by page charge payment. Therefore, and solely to indicate this fact, this article is hereby marked "advertisement" in accordance with 18 USC section 1734.

© 2017 by The American Society of Hematology

CD45, CD38) and the nucleic acid dye Syto41 and gave concordant MRD results with PCR-MRD data in 86.5% of samples. A Chinese study reported an 8-color antibody tube (CD10, CD19, CD20, CD34, CD38, CD45, CD58, plus CD66c or CD13/CD33 or NG2/CD15) with a sensitivity of 0.001% in 81.6% of patients.⁸ Shaver et al elegantly analyzed the relative contribution that each marker and/or pair of markers made to detect MRD¹⁵ and concluded that a single 8-color tube consisting of CD9, CD10, CD19, CD20, CD34, CD38, CD45, and CD58 could provide as much diagnostic utility as their existing 3-tube panel with 12 markers.

Within the EuroFlow Consortium (EU-FP6, LSHB-CT-2006-018708), we aimed to design standardized 8-color immunophenotyping protocols for multicenter MRD measurement in BCP-ALL and to improve the sensitivity of the assay to $\leq 10^{-5}$ (at least comparable to PCR). First, in order to select the most informative markers in distinguishing BCP-ALL from normal BCP cells, we applied novel software tools and principal component-based analyses.^{16,17} In each cycle of design-test-evaluate-redesign, the antibody tubes were tested on BCP-ALL samples and normal and/or regenerating bone marrow (BM), followed by assessment of the contribution of each antibody, until satisfactory results were obtained after 5 testing rounds. Second, a flow cytometric protocol for staining and acquisition of large numbers of cells (>4 million) was developed, allowing theoretical sensitivities of at least 0.001% ($\leq 10^{-5}$). Finally, the selected antibody tubes and standardized laboratory procedures were prospectively validated on follow-up samples from BCP-ALL patients, using the EuroMRD PCR-MRD methods in parallel as gold standard.²

Materials and methods

BCP-ALL patients and normal controls

Data were collected in 7 EuroFlow centers. BM samples obtained from healthy donors or patients in whom no hematological malignancy could be detected (eg, BM samples submitted for lymphoma staging, neuroblastoma staging) were used as control BM for normal/reactive BCP cells. BM samples obtained from pediatric ALL patients after induction therapy (day 78 of therapy) or 1 year after stop of therapy, proven to be MRD-negative by real-time quantitative polymerase chain reaction (RQ-PCR) analysis, were used as a source of regenerating BCP. In the first part of the study (panel design and optimization), samples from 319 BCP-ALL patients, which were consecutively received during 5 design-test-evaluate-redesign phases (initial phase: n = 69; phase 1: n = 61; phase 2: n = 28; phase 3: n = 78; phase 4: n = 83), were included. In the second part of the study (MRD analysis), 377 follow-up samples obtained from 178 BCP-ALL patients (day 15: n = 111; day 33: n = 139; day 78: n = 107; other time points: n = 20) were included. Patient characteristics are summarized in Table 1. The institutional review board of each participating center approved this study, and informed consent for study participation was obtained from each patient and/or his/her legal guardian.

Immunophenotyping MRD panel design

First, BM samples obtained from 69 BCP-ALL patients at diagnosis were stained with the EuroFlow BCP-ALL antibody panel (23 different antibodies in four 8-color tubes).¹⁸ The subsequently designed and optimized MRD tubes were tested during phase 1 to 4 on diagnostic BM samples from BCP-ALL patients using the standardized EuroFlow sample preparation and instrument setup protocols.^{18,19} Data were analyzed using Infinicyt software by comparing BCP-ALL cells with the nearest normal/reactive BCP subsets using automated population separator (APS) plots (see supplemental Methods, available on the *Blood* Web site) as illustrated in Figure 1 and supplemental Figures 1 and 2. Regenerating BCP cells from 6 T-ALL patients were used as an additional negative control (supplemental Figure 3).

Immunophenotyping MRD analyses

The finally selected BCP-ALL MRD tubes were evaluated on BM samples obtained during follow-up of BCP-ALL patients, using an optimized bulk-lysis protocol (see "Results").²⁰ BM samples were processed according to this new EuroFlow bulk-lysis protocol and subsequently stained using the regular EuroFlow protocol.¹⁹ MRD analyses and interpretation were performed locally, and data were subsequently sent to the BCP-ALL-MRD coordinator for central evaluation. Initial FCM-MRD data analysis was performed using 2-dimensional dot plots for sequential gating of BCP-ALL cells, comparable to previous studies using 4 to 6 color stainings.^{10,12} For this study, we provisionally defined a minimum of 10 clustered events to consider a sample as MRD positive (lower limit of detection, LOD) and a minimum of 40 clustered events for accurate quantitation of the MRD level (lower limit of quantitation, LLOQ).²¹ Interlaboratory variability in data analysis was evaluated as described in the supplemental Methods and supplemental Figure 4.

RQ-PCR-based MRD analyses

MRD levels were routinely determined by RQ-PCR analysis of rearranged immunoglobulin and/or T-cell receptor (TR) gene rearrangements in laboratories participating in the quality control rounds of the EuroMRD network (see www.EuroMRD.org).^{3,22-26} RQ-PCR data, performed in triplicate, were analyzed according to the EuroMRD guidelines, using the criteria to prevent false-negative MRD results.² Because application of these criteria might result in some false-positive RQ-PCR results,² we performed next-generation sequencing (NGS) to confirm or exclude the presence of MRD in discordant samples considered positive by RQ-PCR but negative by FCM.

NGS-based MRD analyses

NGS was generally performed as described previously.²⁷ Briefly, depending on the *IGH*, *TRG*, and/or *TRD* rearrangements applied as MRD targets in the RQ-PCR analysis, we performed a targeted approach: the follow-up samples were amplified using the multiplex primer set(s) of the relevant immunoglobulin/TR locus only, and data analysis was focused on the specific junctional region sequence (ie, the 1 used for RQ-PCR analysis). The primers for TCRG were newly designed (supplemental Table 1), and individual primer combinations from multiplex PCR were tested for sensitivity using NGS for diluted diagnostic ALL samples from patients with respective V and Jgamma segment combinations, all reaching the sensitivity of 10^{-5} . All data were finally scored as either MRD-positive or MRD-negative.

Results

Design and optimization of 8-color MRD labeling for BCP-ALL

In the initial phase, 5 antibodies (CD19, CD45, CD34, CD10, and CD20) were upfront selected as backbone markers because they allow appropriate BCP gating as well as characterization of several BCP subpopulations and are known to allow discrimination between normal BCP and BCP-ALL cells.^{10,28-30} To evaluate which other markers could contribute to optimal separation of BCP-ALL cells from normal/reactive BCP cells, the EuroFlow BCP-ALL diagnosis panel¹⁸ was applied to 69 BCP-ALL patients as well as to normal/reactive BM samples. Based on principal component analysis (visualized through APS plots)^{16,17,19} of the BCP-ALL cells vs normal/reactive BCP cells (analyzed per tube), CD9, CD123, CD66c, CD81, CD24, and CD10 appeared to be markers that were most frequently differentially expressed (see supplemental Figure 5). These markers were combined with the 5 backbone markers listed above and complemented with terminal deoxynucleotidyltransferase (TdT) and CD58, both previously reported to be of relevance for BCP-ALL MRD analyses.^{10,28,29,31} The remaining open position was filled in with surface membrane IgKappa/IgLambda (SmIgK/L), as a potential exclusion maker for

Table 1. Patient characteristics

	Initial phase	%	Phase 1	%	Phase 2	%	Phase 3	%	Phase 4	%	MRD phase	%
Number	69		61		28		78		83		178	
Age, y												
Median	6		5		5		5		3		5	
Minimum	0		1		1		0		0		0	
Maximum	76		75		17		17		55		77	
Sex												
Male	35	56	40	66	14	50	46	59	43	52	94	53
Female	27	44	21	34	14	50	32	41	40	48	84	47
White blood cell count (×10⁹/L)												
Median	13.4		8.2		13.7		10.9		15.9		8.5	
Minimum	0.4		0.8		1.2		0.7		0.5		0.5	
Maximum	818		87.6		312		595		1250		900	
No data	10		1		0		2		19		11	
Immunophenotype												
Pro-B-ALL	8	13	2	3	4	14	5	6	6	7	5	3
Common-ALL	44	73	48	79	20	71	59	76	56	67	137	77
Pre-B-ALL	8	13	11	18	4	14	13	17	21	25	35	20
No data	7		0		0		1		0		1	
Genetic data*												
<i>TEL-AML1</i>	6	11	15	27	6	27	11	16	11	18	45	26
<i>BCR-ABL</i>	4	6	1	2	0	0	3	4	1	2	4	2
<i>E2A-PBX1</i>	0	0	3	6	0	0	0	0	0	0	4	2
<i>MLL-AF4</i>	4	6	1	2	2	7	3	4	4	6	3	2
<i>MLL</i> other	2	3	1	2	1	4	1	1	1	2	2	1
Hyperdiploid	14	25	27	48	6	26	26	37	22	37	67	43
Hypodiploid	0	0	1	2	1	4	1	1	1	2	2	1

*Genetic data were not available for all patients; percentages refer to positive patients per all analyzed patients.

more mature BCPs. Fluorochrome positions were primarily determined based on the position of the involved markers in the EuroFlow BCP-ALL panel.¹⁸

The resulting 3 phase 1 MRD tubes (Table 2) were subsequently tested on 61 consecutive BCP-ALL patients at diagnosis, and the discriminatory power was evaluated by comparing the leukemic BCP with the nearest normal BCP subset in APS plots. Whereas both tube 1 and tube 3 gave good/fair separation in ~60% of cases, tube 2 was clearly less informative (fair/good separation in <35% of cases) (Figure 2). When the tube providing the best separation for each patient was selected, good/fair separation was observed in 77% of cases. Considering only tube 1 and 3, good/fair separation was still observed in 71% of cases. These data indicate that tube 1 and 3 had complementary value and confirm the limited value of tube 2.

To evaluate the relevance of each individual marker in discriminating BCP-ALL cells from normal/reactive BCP cells, those markers that received a weight over 10% in the first or second principal component in an APS view of the nearest normal BCP cells and the BCP-ALL cells were selected. CD66c (80% of cases), CD9 (63%), and CD123 (55%) contributed most frequently.

Based on these phase 1 BCP-ALL results, the panel was redesigned: CD58, TdT, SmIgK/L, and CD81 were (at least provisionally) excluded, whereas CD22, which might be important for gating of B-cells in case of CD19-targeting therapies, was included (Table 2).

The phase 2 BCP-ALL MRD tubes were evaluated on diagnostic samples from 28 consecutive BCP-ALL patients. Good/fair separation between BCP-ALL cells and their nearest normal/reactive BCP counterpart was possible in ~75% of cases in both tubes (Figure 2), showing significant improvement over the phase 1 tubes. If the best score of both tubes was used for each case, over 85% of BCP-ALL cases showed good/fair separation from the corresponding normal BCP subset, and in only 3 cases (12%) separation was poor.

During phase 2, additional studies were performed: (1) Because of nonoptimal (relatively weak) CD9 (MEM61) staining, another CD9 clone (ML13) was evaluated with much stronger results; (2) Two newly available fluorochromes (APC-C750 and APC-A750) showed lower background than APC-H7 (less binding to apoptotic cells; no binding to monocytes); (3) More detailed evaluation of the usefulness of CD81 vs CD24 (reevaluation of phase 1 data) showed that CD81 was more frequently differentially expressed between normal/reactive BCP cells and BCP-ALL cells and showed that CD81 in combination with only the backbone markers resulted in a higher percentage of cases with good separation than CD24 did (31% vs 20%); (4) Because CD66c and CD123 are both virtually negative on normal/reactive BCP cells, we tested whether these markers could be combined in the phycoerythrin (PE) channel and concluded that background levels were not affected by combining these 2 markers (data not shown). Based on these data, the new CD9 clone was included in the MRD panel, APC-H7 was replaced by APC-C750/A750, CD24 was replaced by CD81, and CD66c and CD123 were combined into 1 fluorescence channel. The open fluorescein isothiocyanate (FITC) position was used for further evaluation of CD58. The combined data provided the phase-3 BCP-ALL MRD tubes (Table 2).

The 2 phase 3 BCP-ALL MRD tubes were evaluated on 78 BCP-ALL patients. Overall, tube 1 resulted in good/fair separation in 90% of cases, whereas this was achieved in 82% of cases for tube 2 (Figure 2). In the 3 cases for which tube 1 did not result in good separation, normal/reactive BCP and BCP-ALL cells could be separated in tube 2, mainly due to differential expression of CD81. Further evaluation showed that CD38 (~35% of cases), CD66c/CD123 (~30%), and CD81 (~19%) improved the separation between normal/reactive and malignant BCP cells as compared with the 5 backbone markers only, whereas CD9, CD58, and CD22 had no or limited additional value. CD9 in tube 1 was




Figure 1. Data analysis strategy used to optimize the antibody panel for distinguishing between BCP-ALL cells and their nearest normal/reactive BCP counterpart. First, multiple normal/reactive BM samples and/or regenerating BM samples were merged (phase 1: $n = 7$; phase 2: $n = 11$; phase 3: $n = 14$; phase 4: $n = 10$) and CD19-positive B cells were selected. These were subdivided into 4 B-cell subsets, based on the backbone markers (CD19/CD10/CD20/CD34/CD45): CD34⁺ pre-B-I cells (light green), CD34⁻/CD10⁺/CD20⁻ to dim pre-B-II/immature cells (dark blue), CD34⁻/CD10⁺/CD20⁺ immature/transitional B cells (light blue), and CD34⁻/CD10⁻/CD20⁺ mature B cells (dark green). Dot plots of CD34 vs CD10 (A) and CD10 vs CD20 (B) are shown. The 1 standard deviation (SD) (dashed line) and 2 SD lines (solid line) of the 2 most immature BCP subsets (pre-B-I [light green] and pre-B-II/immature [dark blue]) were displayed in an APS view, which was subsequently fixed (supervised; C). Each individual BCP-ALL case was added to the fixed APS plot, and the normal BCP population nearest to each of the BCP-ALL populations was defined (D). The BCP-ALL cells and nearest normal BCP subset were then visualized in a separate (nonfixed and balanced) APS plot, 1 using the backbone markers only (E) and 1 using all 8 markers (F), by plotting the 1 SD curve and 2 SD curves of the 2 populations. To prevent an influence of the number of cells on the principal component analysis (PCA), we opted to use a balanced PCA, implying a fixed ratio between normal and pathological events. Finally, the separation between the 2 populations was scored based on: no overlap between 2 SD curves: 3 points; overlap of the 2 SD curves: 2 points; overlap of the 2 SD and the 1 SD curves: 1 point; overlap of both 1 SD curves: 0 points. An example of this scoring is shown in supplemental Figure 1. It should be noted that the above described strategy was only used for optimizing the antibody panel for the BCP-ALL MRD tubes and not for actual MRD analyses.

therefore replaced by CD81-FITC (which demonstrated equally good staining patterns as CD81-APC-C750) and tube 2 was discarded.

Because in a few cases the evaluated MRD tubes still did not yet result in sufficient separation between normal/reactive and malignant BCP cells, we evaluated several other markers reported to be of potential interest for MRD analysis (eg, CD44-FITC, CD27-PE, CD164-FITC, CD73-PE, CD49f-FITC, CD200-PE, CD86-FITC, and Drebrin-PE).³²⁻³⁶ Based on initial testing on diagnostic BCP-ALL samples, CD73 and CD304 appeared to be most promising based on the level and frequency of overexpression ($\sim 20\%$ for CD73 and $\sim 40\%$ for CD304) and their stability during follow-up (data not shown). Because it appeared not to be possible to combine these 2 markers with CD66c and CD123 in a single fluorescence channel (due to too high background levels), a second tube was designed; this tube was identical to tube 1 but with CD73/CD304 instead of CD66c/CD123 in the PE channel (Table 2).

The 2 phase 4 BCP-ALL MRD tubes were run on 83 consecutive diagnostic BCP-ALL samples. Overexpression of CD66c/CD123 or CD73/CD304 was observed in 45% and 46% of cases, respectively; 31% of cases did not show overexpression of either CD66c/CD123 or

CD73/CD304. Tube 1 resulted in good/fair separation in 89% of cases, whereas this was attained in 82% of cases for tube 2 (Figure 2). If the best score of both tubes was used, 99% of cases showed good/fair separation between the BCP-ALL cells and the nearest normal/reactive BCP subset. Therefore, tube 1 and tube 2 were complementary to each other, and one might either decide at diagnosis which tube is best for monitoring the particular patient or use both tubes to have an extra internal control and more precise measurements. These 2 optimized tubes were considered to be final and ready for further evaluation in follow-up samples of BCP-ALL patients.

Optimization of the flow cytometric MRD sample preparation protocol

We aimed for a sensitivity of $\leq 10^{-5}$, at least comparable to the sensitivity reached in RQ-PCR-based MRD analysis. If a cluster of 10 to 40 BCP-ALL cells should be present to consider a sample as positive, one should acquire at least 4 million cells in order to reach the required sensitivity. Because the cellularity of BM samples obtained during the early phases of treatment is frequently low,³⁷ staining whole BM

Table 2. Development and design of the EuroFlow BCP-ALL MRD panel

Phase Initial*	Violet laser		Blue laser				Red laser	
	PB	PO	FITC	PE	PerCP Cy5.5	PECy7	APC	APC H7
1	CD20	CD45	CD58	CD66c	CD34	CD19	CD10	CD38
	CD9	CD45	TdT	CD13	CD34	CD19	CD22	CD24
	CD21	CD45	CD15/CD65	NG2	CD34	CD19	CD123	CD81
	PB	PO	FITC	PE	PerCP Cy5.5	PECy7	APC	APC H7
	CD20	CD45	CD58	CD66c	CD34	CD19	CD10	CD38
2	CD20	CD45	TdT	SmlgK/L	CD34	CD19	CD10	CD24
	CD20	CD45	CD9	CD123	CD34	CD19	CD10	CD81
	PB	PO	FITC	PE	PerCP Cy5.5	PE Cy7	APC	APC H7
	CD20	CD45	CD9	CD22	CD34	CD19	CD10	CD38
	CD20	CD45	CD123	CD66c	CD34	CD19	CD10	CD24
3	PB	PO	FITC	PE	PerCP Cy5.5	PE Cy7	APC	APC C750
	CD20	CD45	CD9 (ML13)	CD66c/CD123	CD34	CD19	CD10	CD38
	CD20	CD45	CD58	CD22	CD34	CD19	CD10	CD81
4	PB	PO	FITC	PE	PerCP Cy5.5	PE Cy7	APC	APC C750
	CD20	CD45	CD81	CD66c/CD123	CD34	CD19	CD10	CD38
	CD20	CD45	CD81	CD73/CD304	CD34	CD19	CD10	CD38

Markers that were changed as compared with the previous panel are marked in boldface.

*Only tubes 1, 3, and 4 of the EuroFlow BCP-ALL diagnosis panel were evaluated, because most markers in tube 2 (Cy μ , SmlgK, SmlgL, SmlgM) were not expected to contribute to separation of normal and malignant BCP cells.

samples using the regular EuroFlow protocols would not allow acquisition of millions of cells. We therefore designed and tested a new EuroFlow erythrocyte bulk-lysis procedure: sufficiently large volumes of BM, ie, containing >10 million cells, are lysed, and the leukocytes are subsequently resuspended in a small volume of washing buffer. This new protocol allowed staining of 10 million cells in 100 μ L of cell suspension per tube (supplemental Table 2). Evaluation of this new protocol showed that the percentage of doublets did not increase, that the number of evaluable leukocytes increased significantly, and that there were no major differences in cellular composition as compared with the regular EuroFlow staining protocol (Figure 3). Given the large

increase in the number of cells stained with this new approach, all antibody titers were reevaluated; modifications appeared not to be necessary.

Evaluation of the EuroFlow BCP-ALL MRD tube

To evaluate whether the newly designed high-throughput EuroFlow BCP-ALL MRD strategy performed well, we tested the final MRD tubes on follow-up samples of BCP-ALL patients. Based on the immunophenotype of the BCP-ALL cells at diagnosis, 1 MRD tube was selected and subsequently used for MRD evaluation. First, flow




Figure 2. Power to distinguish BCP-ALL cells from their nearest normal BCP counterpart using the EuroFlow BCP-ALL MRD tubes. Data reflect the percentage of patients that reached the specified score, obtained as described in supplemental Figure 1. Briefly, for each patient, BCP-ALL cells and their nearest normal BCP subset were visualized in a (nonfixed and balanced) APS plot showing the median, 1 SD curves, and 2 SD curves for both populations. Each patient was subsequently scored as follows: no overlap between 2 SD curves: 3 points; overlap of the 2 SD curves: 2 points; overlap of the 2 SD and the 1 SD curve: 1 point; overlap of both 1 SD curves: 0 points. Max, the maximal score of the individual tubes. Phase 1: 7 normal/reactive BM samples and 61 BCP-ALL patients; phase 2: 7 normal BM samples, 4 regenerating BM samples, and 28 BCP-ALL patients; phase 3: 5 normal/reactive BM samples, 9 regenerating BM samples, and 78 BCP-ALL patients; phase 4: 10 normal/reactive BM samples and 83 BCP-ALL patients.




Figure 3. Evaluation of the EuroFlow bulk-lysis protocol. For reasons of comparison, each of the BM samples (day 15: $n = 15$; day 33: $n = 15$; day 78: $n = 12$) was processed according to the standard EuroFlow protocol (FL) and in parallel according to the EuroFlow bulk-lysis protocol (BL). (A) Number of leukocytes, debris, and doublets, calculated as percentage of acquired events. Using the bulk-lysis method, significantly less debris ($P = .032$ by paired Student t test) and significantly more leukocytes ($P = .03$ by paired Student t test) were measured. There were no significant differences between the 2 methods for the percentage of doublets. (B) Absolute number of leukocytes acquired. Using BL, on average 12-fold more leukocytes could be acquired ($P < .0001$). Please note that we included relatively many day 15 samples in order to be able to evaluate the impact of the 2 methods on the MRD levels as well. However, these day 15 samples generally have a very low white blood cell count; consequently, the number of leukocytes acquired after BL is still relatively low in a subset of samples. (C) Distribution of leukocyte subpopulations, defined as percentage of leukocytes. By paired Student t test (2-sided), small but statistically significant differences were observed for T/NK cells (mean: 24% vs 26%, $P = .0047$), granulocytes (mean: 33% vs 38%, $P < .001$), and monocytes (mean: 3.2% vs 4.5%, $P < .001$), whereas no significant differences were observed for the remaining populations. Of note, in 2 samples, MRD was only detected using the bulk-lysis method (0.013% and 0.018%) but not using the whole BM method. In the 11 samples MRD positive by both methods, MRD levels were not significantly different from each other (paired Student t test: $P = .30$), with mean values of 6.3% and 6.7% by whole BM and bulk-lysed BM method, respectively. Correlation analysis showed a Spearman r of 0.964 (95% confidence interval: 0.857-0.991; $P < .0001$). PC, plasma cells.

cytometric MRD data obtained in 178 BCP-ALL patients were compared with routinely obtained PCR-MRD data. As shown in Figure 4A, the concordance between the FCM-MRD data and PCR-based MRD data was highly dependent on the number of cells acquired by FCM. In addition, the sensitivity of FCM-MRD (percentage samples positive by both FCM and PCR relative to samples positive by PCR) significantly increased when higher cell numbers were acquired (Figure 4B). Therefore, only samples in which MRD could clearly be detected by FCM-MRD or samples that had sufficient cells acquired for reaching a sensitivity of $\leq 10^{-5}$ were included in the subsequent analyses. Based on an LLOQ of 40 events, at least 4×10^6 cells should be acquired, which was possible in 227 out of 377 samples (60%). FCM-MRD data obtained in these patients were comparable to PCR-based MRD results in 93% of samples (Figure 5; Table 3). All but 1 of the 17 discordant samples (7 FCM⁺/PCR⁻ and 10 FCM⁻/PCR⁺) had MRD levels $< 10^{-4}$ (supplemental Table 2). Bland-Altman analysis showed higher PCR-based MRD values with a mean difference of 0.34 log or a factor of 2.2 (supplemental Figure 6).

Detailed evaluation of discordant cases

To evaluate the discordant cases, several additional analyses were performed. First, FCM fetal calf serum data files were blindly distributed to 4 laboratories for reanalysis of the FCM-MRD data. Out of the 7 cases initially scored positive by FCM-MRD and negative by PCR-MRD, 6 were interpreted as negative by all 4 centers upon FCM-MRD reanalysis, whereas 1 sample (day 15) was consistently scored positive by all 4 centers. Second, RQ-PCR-MRD data were checked for cases negative by FCM-MRD and positive by PCR-MRD.

In 8 of 10 cases (all confirmed to be FCM-MRD–negative by reanalysis in different centers), PCR-MRD data were considered positive based on a single well in a single target. To further evaluate whether this low-level positivity was potentially caused by nonspecific amplification, NGS was used to confirm the possible presence of the leukemia-specific antigen receptor rearrangement. In 7 out of 9 available samples, NGS-MRD was negative, whereas MRD positivity could be confirmed in the remaining 2 patients. Thus, 6 of 17 (35%) discordant cases were due to initial misinterpretation of the FCM-MRD data and at least 7 of 17 (41%) were due to overinterpretation of PCR-MRD data; the remaining 4 cases appeared to be truly discordant cases (supplemental Table 3). After these additional evaluations, the actual concordance increased to 98%. If only samples with MRD levels $< 0.01\%$ were included, 97% gave concordant results.

Discussion

After 5 phases of optimization, we finally selected two 8-color antibody tubes that only differed for the markers present in the PE channel (CD66c/CD123 vs CD73/CD304). These 2 tubes are comparable with 8-color BCP-ALL MRD tubes recently used in other studies, because all proposed panels include CD19, CD10, CD20, CD34, and CD45. In our study, these markers were considered backbone markers from start onwards, based on our previous experience.^{10,12,18,28-30} Also, CD38 is present in all proposed panels and was proven to be relevant in our present study as well as in the studies of Karawajew et al¹⁴ and Shaver et al.¹⁵ The remaining 2 positions were completed with different




Figure 4. Performance of FCM-MRD vs PCR-based MRD is dependent on the number of acquired cells. (A) The percentage of discordant cases by FCM-MRD and PCR-MRD is shown for individual time points (day 15, day 33, and day 78) as well as for all samples together. Data are presented for variable numbers of acquired cells: all samples (independent of cell number) and samples with at least 1, 2, 3, 4, or 5 million cells acquired. (B) The sensitivity of FCM-MRD relative to PCR-MRD is shown for individual time points (day 15, day 33, and day 78) as well as for all samples together. Data are presented for variable numbers of acquired cells: all samples (independent of cell number; $n = 377$) and samples with at least 1 ($n = 330$), 2 ($n = 287$), 3 ($n = 255$), 4 ($n = 227$), or 5 million cells ($n = 191$) acquired. Sensitivity is calculated as the number of samples positive by both FCM and PCR divided by the total number of samples positive by PCR (ie, the reference method).

markers: CD9, CD13/CD33, CD15/NG2, CD58 (2 studies), CD66c (2 studies), CD73, CD81, CD123, CD304, and a nucleic acid dye. In our analysis, CD9, CD58, and CD22 appeared to be of limited value and therefore were discarded, while Shaver and colleagues, who also applied mathematical modeling systems, identified these as important MRD markers.¹⁵ However, they did not test CD66c, CD73, CD123, or CD304, and the difference between the contribution of CD9 and CD81 in their study was limited.¹⁵ CD15/NG2 might be relevant in ALL with MLL gene rearrangements, mainly occurring in infants.³⁸ These cases are rare and frequently present with a pro-B-ALL immunophenotype that can relatively easily be distinguished from normal BCP cells and plasma cells.³⁹ In our study, we tested and finally selected 4 markers that are frequently abnormally expressed on BCP-ALL cells: CD66c (associated with BCR-ABL and hyperdiploidy),^{40,41} CD123 (associated with hyperdiploidy⁴²), CD73,³⁶ and CD304 (possibly associated with TEL-AML1).³³ By combining 2 of these markers in a single fluorescence channel, abnormal expression could be identified in ~70% of BCP-ALL patients. Furthermore, in combination with the backbone markers and CD81, BCP-ALL cells could clearly be distinguished from normal BCP cells in virtually all patients. Thus, after multiple phases of multicenter testing of a wide range of leukocyte markers, antibody clones, and fluorochrome-conjugated reagents using objective novel software tools, we were able to select 2 highly effective BCP-ALL MRD tubes.

It remains to be evaluated whether the designed BCP-ALL MRD tubes can also be used during antibody-based therapies. Especially Blinatumomab and chimeric antigen receptor-T-cells (targeting CD19) may hamper the gating of BCP based on CD19. Although alternative gating strategies can be applied (eg, based on CD10, CD34, and/or CD45), one could also decide to add CD24 and/or CD22 to the current tubes (transforming it into a 10-color tube) in MRD-based trials involving Blinatumomab.⁴³ Addition of CD24 and CD22 will also have the advantage that the earliest BCP cells, expressing CD24 and/or

CD22 but not yet CD19,⁴⁴ can be identified; this may be of relevance for the identification of all BCP cells in regenerating BM samples.

In order to obtain MRD data with good sensitivity, acquisition of large numbers of cells appears to be a prerequisite. There is no consensus yet about the number of cells needed for a population. Most studies in BCP-ALL indicate a minimum number between 10 and 50 events, whereas a recent consensus report on MRD detection in multiple myeloma patients defined 20 and 50 cells as the LOD and LLOQ, respectively.^{14,21,45} Consequently, a sensitivity of 10^{-5} (generally reached in PCR-MRD and NGS-MRD analysis) requires acquisition of $\geq 10^6$ cells, preferably $\geq 5 \times 10^6$ cells. We therefore developed the new Euroflow bulk lysis protocol, allowing acquisition of such high cell numbers. Although the bulk lysis protocol contains several washing steps, which will likely result in some cell loss, there is no evidence for selective loss of BCP-ALL cells, given the high concordance between the final FCM-MRD results and the PCR-MRD results. To our best knowledge, other FCM-MRD studies so far have not acquired $\geq 4 \times 10^6$ events and therefore could not have reached the same sensitivity as shown here, although the study by the ALL-REZ-BFM 2002 trial group comes close.¹⁴ Our data clearly show that acquisition of large numbers of cells ($\geq 4 \times 10^6$) is a prerequisite for obtaining good sensitivities and data that are truly comparable to PCR-MRD data.

We finally tested the newly designed and optimized BCP-ALL MRD tubes in combination with the bulk lysis protocol on follow-up samples from BCP-ALL patients and compared the FCM-MRD data with PCR-MRD data. The concordance between both methods (in the absence of any cutoff) was extremely high (93%) and was significantly better than in previous studies (82.3%¹⁰ and 86.5%¹⁴). As mentioned above, this increased concordance is likely due to the higher sensitivity of the current study, resulting from the higher number of cells analyzed. Detailed evaluation of the discordant samples using NGS-based approaches showed that several discrepant cases were




Figure 5. Comparison of MRD data obtained by 8-color EuroFlow FCM and routinely obtained molecular MRD data. Flow cytometric MRD data were compared with molecular MRD data and are shown for samples obtained at day 15 (A), day 33 (B), or day 78 (C). In the lower right part of each panel, the number of FCM⁺/PCR⁺, FCM⁺/PCR⁻, FCM⁻/PCR⁺, and FCM⁻/PCR⁻ is indicated. Only samples in which MRD could clearly be detected by FCM-MRD or samples which had sufficient cells acquired for reaching a sensitivity of $\leq 10^{-5}$ (ie, $\geq 4 \times 10^6$ cells acquired) were included in the analyses. Based on a LOD of 10 events, the sensitivity of FCM-MRD was 2.5×10^{-6} ; the quantitative range (based on a LLOQ of 40 events) is 10^{-5} . Consequently, FCM-MRD data between 2.5×10^{-6} and 10^{-5} should be considered positive, but are below the limit of quantitation.

Table 3. Concordance between flow cytometric and molecular MRD data

	Day 15	Day 33	Day 78	Other TP	Total	% (Total %)
Concordant						
FCM ⁺ /PCR ⁺	102	47	11	9	169	74 (93)
FCM ⁻ /PCR ⁻	1	5	31	4	41	18
Discordant						
FCM ⁻ /PCR ⁺	0	5	5	0	10	4 (7)
FCM ⁺ /PCR ⁻	1	4	2	0	7	3
Total	104	61	49	13	227	100

Only samples in which MRD could clearly be detected by FCM-MRD or samples which had at least 4×10^6 cells acquired were included in the analyses. Based on an LLOQ of 40 events, this allowed a sensitivity of at least 10^{-5} .

TP, time points.

due to overinterpretation of the PCR-MRD data, as the involved leukemia-specific immunoglobulin/TR rearrangements could not be detected by (qualitative) NGS analysis. Consequently, it is most likely that in these discordant cases the positive RQ-PCR MRD data, interpreted according to the EuroMRD criteria for prevention of false-negative MRD results, are due to nonspecific amplification,^{2,10,46,47} and that these samples actually were MRD-negative.

Next to these “false-positive PCR-MRD” samples, part of the initially discordant cases could be explained by “false-positive FCM-MRD” results: samples were initially scored MRD-positive (generally at very low levels of $<0.01\%$), but these samples were consistently scored MRD-negative upon blind reanalysis at 4 different centers. Interpretation of FCM-MRD data, especially at MRD levels $<0.01\%$, is still expert-based and depends on the number of events in the suspected population, their distance from normal, and the homogeneity of the suspected population (clustering of suspected cells). Although the number of events may easily be defined, distance from normal and homogeneity of the population are more complex to be objectively defined. Novel software tools, including automated gating approaches and maturation pathway analysis (Supplemental Figure 7) are currently being developed within the EuroFlow Consortium and will facilitate more standardized and objective FCM-MRD measurements in the near future.³⁷ Preliminary maturation-based FCM-MRD data showed very good concordance with PCR-MRD data, although further improvements are needed for detection of low levels of MRD ($<0.01\%$).

The NGS-MRD analyses and reanalyses of FCM-MRD data increased the concordance to an unprecedented rate of 98%; the remaining discordant cases are likely due to statistical variation around the detection limits of both assays.¹⁴ Therefore, the here presented EuroFlow FCM-MRD strategy proved to be highly sensitive (at least comparable to PCR-MRD) and fast and allows standardized quantification of MRD in virtually all BCP-ALL patients. By increasing the number of acquired cells to 10^7 (eg, by running both tubes with 5×10^6 cells), the sensitivity and robustness likely can even be further increased.

Acknowledgments

This research was performed within the EuroFlow Consortium, which started with an EU-FP6 grant (LSHB-CT-2006-018708) and obtained sustainability by protecting and licensing intellectual property, thereby obtaining royalties, which are exclusively being used for supporting the EuroFlow research program (chairmen: J.J.M.v.D. and A.O.). The authors gratefully acknowledge Tomas Kalina and the technicians of the Laboratory for Medical Immunology, Erasmus MC, for technical

assistance and thank Christa Homburg and colleagues for performing molecular MRD analyses of part of the Dutch patients and Marieke Comans-Bitter for organizational support.

E.F. was supported by the Grant Agency of the Czech Republic (project of Centre of Excellence No. P302/12/G101). L.S., T.S., and P.T. were supported by European Research Area Network (ERA-NET) PrioMedChild, grant 40-41800-98-027. E.M. was supported by Ministry of Health of the Czech Republic, grant 15-28525A. M.K. was supported by the University Hospital Motol, Prague, Czech Republic (00064203). E.S.d.C., Q.L., and A.O. acknowledge the Bilateral Cooperation Program between Coordenação de Aperfeiçoamento de Pessoal de Nível Superior (Brasília, Brazil) and Dirección General de Políticas Universitarias–Ministerio de Educación, Cultura y Deportes (Madrid, Spain) (311/15). E.S.d.C. acknowledges Research Foundation of the State of Rio de Janeiro, Rio de Janeiro, Brazil (E26/110.105/2014, E26/102.191/2013) and Conselho Nacional de Desenvolvimento Científico e Tecnológico–CNPQ of Brazil (400194/2014-7). G.G., C.B., and P.B. were supported by Fondazione Tettamanti. The authors gratefully acknowledge the contribution of the EuroClonality/EuroMRD NGS network (chair: A. W. Langerak) for support with the NGS data. The research for this manuscript was in part performed within the framework of the Erasmus Postgraduate School Molecular Medicine.

Authorship

Contribution: V.H.J.v.d.V., E.M., J.J.M.v.D., and A.O. designed the research; P.T., E.M., Q.L., G.E.G., A.O., J.J.M.v.D., and V.H.J.v.d.V. developed the methodology and the research strategies; P.T. and V.H.J.v.d.V. interpreted the data; P.T., E.M., L.S., A.J.v.d.S.-G., C.B., E.O., G.G., M. Bartels, E.S.d.C., M.N., P.B., T.S., L.L., G.E.G., A.O., and V.H.J.v.d.V. analyzed the FCM data; P.T., E.M., L.S., A.J.v.d.S.-G., M. Bartels, M.N., E.S., O.H., E.S.d.C., J.G.t.M., C.B., E.O., P.B., and L.L. performed FCM experiments; M.K., E.F., J.T., and M. Brüggemann designed and performed the NGS experiments and analyzed NGS data; and V.H.J.v.d.V. wrote the paper; all authors reviewed and approved the manuscript.

Conflict-of-interest disclosure: J.J.M.v.D. has performed contract research for Roche, Amgen, and BD Biosciences; V.H.J.v.d.V. has performed contract research for Roche, Amgen, Pfizer, Janssen, and BD Biosciences and received consultancy fees from Celgene. G.E.G. is an employee of Cytognos SL, Salamanca, Spain. The remaining authors declare no competing financial interests.

A complete list of the members of the EuroFlow Consortium appears in “Appendix.”

Correspondence: Jacques J. M. van Dongen, Department of Immunology, Erasmus MC, University Medical Center Rotterdam, Wytemaweg 80, 3015 CN Rotterdam, The Netherlands; e-mail: j.j.m.van_dongen@lumc.nl.

Appendix: study group members

The members of the EuroFlow Consortium are: Coordination: J. J. M. van Dongen, W. M. Bitter. Participants: V. H. J. van der Velden, A. W. Langerak, P. Theunissen, J. te Marvelde, J. Schilperoord-Vermeulen, A. E. Bras, M. van der Burg, M. Wentink, S. Posthumus (Erasmus MC, Rotterdam, The Netherlands), A. Orfao, J. Almeida, Q. Lecrevisse, J. Flores-Montero,

S. Matarraz, María-Belén Vidriales, J. J. Pérez-Morán, N. Puig, M. Pérez-Andrés, L. Martín, E. Blanco (USAL, Salamanca, Spain), M. Gomes da Silva, A. Medina Almeida, J. Caetano, T. Faria (IPO, Lisbon, Portugal), M. Kneba, M. Brüggemann, M. Ritgen, S. Böttcher, M. Szczepanowski, E. Harbst, M. Ipsen (UNIKIEL, Kiel, Germany), E. Macintyre, V. Asnafi, A. Trinquand, L. Lhermitte (AP-HP, Paris, France), J. Trka, O. Hrusak, T. Kalina, E. Mejstrikova, D. Thurner, M. Novakova (DPH/O, Prague, Czech Republic), T. Szczepanski, L. Sędek, J. Bulsa, A. Sonsala (SUM, Zabrze, Poland), E. Sobral da Costa, C. E. Pedreira, F. Pinto Mariz, T. Sigaud S. Palmereira (UFRJ, Rio de Janeiro, Brazil), M. Lima, A. H. Santos (CHP, Porto, Portugal), E. Sonneveld,

A. J. van der Sluijs-Gelling (DCOG, The Hague, The Netherlands). Affiliated participants: G. Gaipa, C. Buracchi, C. Bugarin (HSGerardo, Monza, Italy), M. Roussel (CHUP, Rennes, France), L. Campos-Guyotat, C. Aanei (CHU, Saint-Etienne, France), N. Villamor (UB, Barcelona, Spain), J. Philippé, C. Bonroy, B. Denys, A. Willems, S. de Sy (UZG, Ghent, Belgium), P. Fernandez (KSA, Aarau, Switzerland), M. Vlkova, J. Nechvatalova (FNUSA, Brno, Czech Republic), A. E. Sousa, A. L. Serra-Caetano (IMM, Lisboa, Portugal), O. Speer, J. Trück (KISPI, Zürich, Switzerland). Associated SME: CYTOGNOS, Salamanca, Spain (M. Martín-Ayuso, J. Bensadon, F. Martin de Lara, G. Grigore, P. Peñalosa, J. Verde, R. Fluxa, P. Manzano).

References

- Flohr T, Schrauder A, Cazzaniga G, et al; International BFM Study Group (I-BFM-SG). Minimal residual disease-directed risk stratification using real-time quantitative PCR analysis of immunoglobulin and T-cell receptor gene rearrangements in the international multicenter trial AIEOP-BFM ALL 2000 for childhood acute lymphoblastic leukemia. *Leukemia*. 2008;22(4):771-782.
- van der Velden VH, Cazzaniga G, Schrauder A, et al; European Study Group on MRD detection in ALL (ESG-MRD-ALL). Analysis of minimal residual disease by Ig/TCR gene rearrangements: guidelines for interpretation of real-time quantitative PCR data. *Leukemia*. 2007;21(4):604-611.
- van der Velden VH, van Dongen JJ. MRD detection in acute lymphoblastic leukemia patients using Ig/TCR gene rearrangements as targets for real-time quantitative PCR. *Methods Mol Biol*. 2009;538:115-150.
- Basso G, Veltroni M, Valsecchi MG, et al. Risk of relapse of childhood acute lymphoblastic leukemia is predicted by flow cytometric measurement of residual disease on day 15 bone marrow. *J Clin Oncol*. 2009;27(31):5168-5174.
- Borowitz MJ, Wood BL, Devidas M, et al. Prognostic significance of minimal residual disease in high risk B-ALL: a report from Children's Oncology Group study AALL0232. *Blood*. 2015;126(8):964-971.
- Cheng SH, Lau KM, Li CK, et al. Minimal residual disease-based risk stratification in Chinese childhood acute lymphoblastic leukemia by flow cytometry and plasma DNA quantitative polymerase chain reaction. *PLoS One*. 2013;8(7):e69467.
- Stow P, Key L, Chen X, et al. Clinical significance of low levels of minimal residual disease at the end of remission induction therapy in childhood acute lymphoblastic leukemia. *Blood*. 2010;115(23):4657-4663.
- Weng XQ, Shen Y, Sheng Y, et al. Prognostic significance of monitoring leukemia-associated immunophenotypes by eight-color flow cytometry in adult B-acute lymphoblastic leukemia. *Blood Cancer J*. 2013;3:e133.
- Dworzak MN, Fröschl G, Printz D, et al; Austrian Berlin-Frankfurt-Münster Study Group. Prognostic significance and modalities of flow cytometric minimal residual disease detection in childhood acute lymphoblastic leukemia. *Blood*. 2002;99(6):1952-1958.
- Denys B, van der Sluijs-Gelling AJ, Homburg C, et al. Improved flow cytometric detection of minimal residual disease in childhood acute lymphoblastic leukemia. *Leukemia*. 2013;27(3):635-641.
- Ryan J, Quinn F, Meunier A, et al. Minimal residual disease detection in childhood acute lymphoblastic leukaemia patients at multiple time-points reveals high levels of concordance between molecular and immunophenotypic approaches. *Br J Haematol*. 2009;144(1):107-115.
- Gaipa G, Cazzaniga G, Valsecchi MG, et al. Time point-dependent concordance of flow cytometry and real-time quantitative polymerase chain reaction for minimal residual disease detection in childhood acute lymphoblastic leukemia. *Haematologica*. 2012;97:1582-1593.
- Thörn I, Forestier E, Botling J, et al. Minimal residual disease assessment in childhood acute lymphoblastic leukaemia: a Swedish multi-centre study comparing real-time polymerase chain reaction and multicolour flow cytometry. *Br J Haematol*. 2011;152(6):743-753.
- Karawajew L, Dworzak M, Ratei R, et al. Minimal residual disease analysis by eight-color flow cytometry in relapsed childhood acute lymphoblastic leukemia. *Haematologica*. 2015;100(7):935-944.
- Shaver AC, Greig BW, Mosse CA, Seegmiller AC. B-ALL minimal residual disease flow cytometry: an application of a novel method for optimization of a single-tube model. *Am J Clin Pathol*. 2015;143(5):716-724.
- Pedreira CE, Costa ES, Almeida J, et al; EuroFlow Consortium. A probabilistic approach for the evaluation of minimal residual disease by multiparameter flow cytometry in leukemic B-cell chronic lymphoproliferative disorders. *Cytometry A*. 2008;73A(12):1141-1150.
- Pedreira CE, Costa ES, Barrena S, et al; EuroFlow Consortium. Generation of flow cytometry data files with a potentially infinite number of dimensions. *Cytometry A*. 2008;73(9):834-846.
- van Dongen JJ, Lhermitte L, Böttcher S, et al; EuroFlow Consortium (EU-FP6, LSHB-CT-2006-018708). EuroFlow antibody panels for standardized n-dimensional flow cytometric immunophenotyping of normal, reactive and malignant leukocytes. *Leukemia*. 2012;26(9):1908-1975.
- Kalina T, Flores-Montero J, van der Velden VH, et al; EuroFlow Consortium (EU-FP6, LSHB-CT-2006-018708). EuroFlow standardization of flow cytometer instrument settings and immunophenotyping protocols. *Leukemia*. 2012;26(9):1986-2010.
- Aalbers AM, van den Heuvel-Eibrink MM, Baumann I, et al. Bone marrow immunophenotyping by flow cytometry in refractory cytopenia of childhood. *Haematologica*. 2015;100(3):315-323.
- Arroz M, Came N, Lin P, et al. Consensus guidelines on plasma cell myeloma minimal residual disease analysis and reporting. *Cytometry B Clin Cytom*. 2016;90(1):31-39.
- van der Velden VHJ, Willemse MJ, van der Schoot CE, Hählen K, van Wering ER, van Dongen JJM. Immunoglobulin kappa deleting element rearrangements in precursor-B acute lymphoblastic leukemia are stable targets for detection of minimal residual disease by real-time quantitative PCR. *Leukemia*. 2002;16(5):928-936.
- van der Velden VHJ, Wijkhuijs JM, Jacobs DC, van Wering ER, van Dongen JJM. T cell receptor gamma gene rearrangements as targets for detection of minimal residual disease in acute lymphoblastic leukemia by real-time quantitative PCR analysis. *Leukemia*. 2002;16(7):1372-1380.
- van der Velden VH, de Bie M, van Wering ER, van Dongen JJ. Immunoglobulin light chain gene rearrangements in precursor-B-acute lymphoblastic leukemia: characteristics and applicability for the detection of minimal residual disease. *Haematologica*. 2006;91(5):679-682.
- Szczepanski T, van der Velden VHJ, Hoogeveen PG, et al. Vdelta2-Jalpha rearrangements are frequent in precursor-B-acute lymphoblastic leukemia but rare in normal lymphoid cells. *Blood*. 2004;103(10):3798-3804.
- Brüggemann M, van der Velden VHJ, Raff T, et al. Rearranged T-cell receptor beta genes represent powerful targets for quantification of minimal residual disease in childhood and adult T-cell acute lymphoblastic leukemia. *Leukemia*. 2004;18(4):709-719.
- Kotrova M, Muzikova K, Mejstrikova E, et al. The predictive strength of next-generation sequencing MRD detection for relapse compared with current methods in childhood ALL. *Blood*. 2015;126(8):1045-1047.
- Lucio P, Gaipa G, van Lochem EG, et al; BIOMED-I. BIOMED-I concerted action report: flow cytometric immunophenotyping of precursor B-ALL with standardized triple-stainings. BIOMED-I Concerted action investigation of minimal residual disease in acute leukemia: international standardization and clinical evaluation. *Leukemia*. 2001;15(8):1185-1192.
- Lúcio P, Parreira A, van den Beemd MW, et al. Flow cytometric analysis of normal B cell differentiation: a frame of reference for the detection of minimal residual disease in precursor-B-ALL. *Leukemia*. 1999;13(3):419-427.
- Mejstriková E, Fronková E, Kalina T, et al. Detection of residual B precursor lymphoblastic leukemia by uniform gating flow cytometry. *Pediatr Blood Cancer*. 2010;54(1):62-70.
- Veltroni M, De Zen L, Sanzari MC, et al; I-BFM-ALL-FCM-MRD-Study Group. Expression of CD58 in normal, regenerating and leukemic bone marrow B cells: implications for the detection of minimal residual disease in acute lymphocytic leukemia. *Haematologica*. 2003;88(11):1245-1252.
- Coustan-Smith E, Song G, Clark C, et al. New markers for minimal residual disease detection in acute lymphoblastic leukemia. *Blood*. 2011;117(23):6267-6276.
- Solly F, Angelot F, Garand R, et al. CD304 is preferentially expressed on a subset of B-lineage

- acute lymphoblastic leukemia and represents a novel marker for minimal residual disease detection by flow cytometry. *Cytometry A*. 2012; 81(1):17-24.
34. Vaskova M, Kovac M, Volna P, et al. High expression of cytoskeletal protein drebrin in TEL/AML1pos B-cell precursor acute lymphoblastic leukemia identified by a novel monoclonal antibody. *Leuk Res*. 2011;35(8):1111-1113.
 35. Vaskova M, Mejstrikova E, Kalina T, et al. Transfer of genomics information to flow cytometry: expression of CD27 and CD44 discriminates subtypes of acute lymphoblastic leukemia. *Leukemia*. 2005;19(5):876-878.
 36. Wang W, Gao L, Li Y, et al. The application of CD73 in minimal residual disease monitoring using flow cytometry in B-cell acute lymphoblastic leukemia. *Leuk Lymphoma*. 2016;57(5):1174-1181.
 37. van Dongen JJ, van der Velden VH, Brüggemann M, Orfao A. Minimal residual disease diagnostics in acute lymphoblastic leukemia: need for sensitive, fast, and standardized technologies. *Blood*. 2015;125(26):3996-4009.
 38. Behm FG, Smith FO, Raimondi SC, Pui CH, Bernstein ID. Human homologue of the rat chondroitin sulfate proteoglycan, NG2, detected by monoclonal antibody 7.1, identifies childhood acute lymphoblastic leukemias with t(4;11)(q21;q23) or t(11;19)(q23;p13) and MLL gene rearrangements. *Blood*. 1996;87(3):1134-1139.
 39. Jansen MW, Corral L, van der Velden VH, et al. Immunobiological diversity in infant acute lymphoblastic leukemia is related to the occurrence and type of MLL gene rearrangement. *Leukemia*. 2007;21(4):633-641.
 40. Kalina T, Vaskova M, Mejstrikova E, et al. Myeloid antigens in childhood lymphoblastic leukemia: clinical data point to regulation of CD66c distinct from other myeloid antigens. *BMC Cancer*. 2005;5:38-49.
 41. Kiyokawa N, Iijima K, Tomita O, et al. Significance of CD66c expression in childhood acute lymphoblastic leukemia. *Leuk Res*. 2014;38(1):42-48.
 42. Djokic M, Björklund E, Blennow E, Mazur J, Söderhäll S, Porwit A. Overexpression of CD123 correlates with the hyperdiploid genotype in acute lymphoblastic leukemia. *Haematologica*. 2009; 94(7):1016-1019.
 43. Cherian S, Miller V, McCullouch V, Dougherty K, Fromm JR, Wood BL. A novel flow cytometric assay for detection of residual disease in patients with B-lymphoblastic leukemia/lymphoma post anti-CD19 therapy. *Cytometry B Clin Cytom*. In press.
 44. van Lochem EG, van der Velden VHJ, Wind HK, te Marvelde JG, Westerdaal NAC, van Dongen JJM. Immunophenotypic differentiation patterns of normal hematopoiesis in human bone marrow: reference patterns for age-related changes and disease-induced shifts. *Cytometry B Clin Cytom*. 2004;60(1):1-13.
 45. Wood BL. Principles of minimal residual disease detection for hematopoietic neoplasms by flow cytometry. *Cytometry B Clin Cytom*. 2016;90(1): 47-53.
 46. Fronkova E, Muzikova K, Mejstrikova E, et al. B-cell reconstitution after allogeneic SCT impairs minimal residual disease monitoring in children with ALL. *Bone Marrow Transplant*. 2008;42(3):187-196.
 47. van der Velden VH, Wijkhuijs JM, van Dongen JJ. Non-specific amplification of patient-specific Ig/TCR gene rearrangements depends on the time point during therapy: implications for minimal residual disease monitoring. *Leukemia*. 2008; 22(3):641-644.