

risk significantly higher than that in matched healthy subjects. Our findings suggest there may be value in regular (perhaps yearly) interval complete blood count follow-up of these patients, analogous to monoclonal gammopathy of undetermined significance (given comparable risks of progression). In addition, consideration of more systematic upfront workup, such as peripheral blood molecular testing for inherited and acquired myeloid mutations, may be of value. Directed study would certainly be required to evaluate and to confirm the optimal initial evaluation and longitudinal surveillance strategy in this patient population.

## Acknowledgments

H.A.-S. is the recipient of the American Society of Hematology Scholar Award.

This work will be included in a dissertation submitted to the Faculty of Harvard Medical School in partial fulfillment of the requirements for the Degree of Master of Medical Sciences in Clinical Investigation (MMSCI) for N.A.

## Authorship

Contribution: N.A. contributed to study design, data collection, data analysis, creation of tables and figures, writing of the first draft of the manuscript, critical revision of the manuscript, and final approval; R.F.G. and D.J.K. contributed to critical revision of the manuscript and final approval; and H.A.-S. contributed to study design, data collection, data analysis, creation of tables and figures, critical revision of the manuscript, and final approval.

Conflict-of-interest disclosure: None of the authors have disclosures relevant to the content of the manuscript. Universal disclosures include the following. D.J.K. has received research funding from Argenx, Biocryst, Immunovant, Rigel, Sanofi (Principia), Takeda (Bioverativ), and UCB; has been a consultant for Alexion (Syntimmune), Amgen, argenx, BioCryst, Bristol Myers Squibb, Caremark, Cellularity, Cellphire, CRICO, Daiichi Sankyo, Hengrui, Immunovant, Incyte, Kyowa-Kirin, Merck Sharp Dohme, Momenta, Novartis, Pfizer, Platelet Biogenesis, Platelet Disorder Support Association, Rigel, Sanofi (Bioverativ), Sanofi (Principia), Sanofi (Genzyme), Sobi (Dova), Takeda, UCB, and Up-To-Date; and owns stock in Rubius. R.F.G. has received research funding from Agios, Dova/Sobi, and Novartis; and has been a consultant for Agios and Principia/Sanofi. H.A.-S. has received research funding from Agios, Dova/Sobi, and Amgen; and been a consultant for Agios, Dova/Sobi, Forma, argenx, Rigel, Moderna, and Novartis. N.A. declares no competing financial interests.

ORCID profiles: N.A., 0000-0002-1482-7792; R.F.G., 0000-0001-7302-0449; H.A.-S., 0000-0001-6175-1383.

Correspondence: Hanny Al-Samkari, Division of Hematology Oncology, Massachusetts General Hospital, Zero Emerson Pl, Suite 118, Office 112, Boston, MA 02114; email: hal-samkari@mgh.harvard.edu.

## Footnotes

Submitted 21 April 2022; accepted 16 August 2022; prepublished online on *Blood* First Edition 26 August 2022.

The online version of this article contains a data supplement.

## REFERENCES

1. Neunert C, Terrell DR, Arnold DM, et al. American Society of Hematology 2019 guidelines for immune thrombocytopenia. *Blood Adv*. 2019;3(23):3829-3866.
2. Ruggeri M, Fortuna S, Rodeghiero F. Heterogeneity of terminology and clinical definitions in adult idiopathic thrombocytopenic purpura: a critical appraisal from a systematic review of the literature. *Haematologica*. 2008;93(1):98-103.
3. Provan D, Arnold DM, Bussel JB, et al. Updated international consensus report on the investigation and management of primary immune thrombocytopenia. *Blood Adv*. 2019;3(22):3780-3817.
4. Rodeghiero F, Stasi R, Gernsheimer T, et al. Standardization of terminology, definitions and outcome criteria in immune thrombocytopenic purpura of adults and children: report from an international working group. *Blood*. 2009;113(11):2386-2393.
5. Adibi P, Faghih Imani E, Talaei M, Ghanei M. Population-based platelet reference values for an Iranian population. *Int J Lab Hematol*. 2007;29(3):195-199.
6. Bain BJ. Ethnic and sex differences in the total and differential white cell count and platelet count. *J Clin Pathol*. 1996;49(8):664-666.
7. Izzi B, Bonaccio M, de Gaetano G, Cerletti C. Learning by counting blood platelets in population studies: survey and perspective a long way after Bizzozero. *J Thromb Haemost*. 2018;16(9):1711-1721.
8. Burrows RF, Kelton JG. Incidentally detected thrombocytopenia in healthy mothers and their infants. *N Engl J Med*. 1988;319(3):142-145.
9. Stasi R, Amadori S, Osborn J, Newland AC, Provan D. Long-term outcome of otherwise healthy individuals with incidentally discovered borderline thrombocytopenia. *PLoS Med*. 2006;3(3):e24.

<https://doi.org/10.1182/blood.2022016757>

© 2022 by The American Society of Hematology

## TO THE EDITOR:

# Improved semiautomated detection of TRBC-restricted Sézary cells unveils a spectrum of clonal cluster immunophenotypes

Pedro Horna, Gregory Otteson, Min Shi, Jansen N. Seheult, Dragan Jevremovic, and Horatiu Olteanu

Division of Hematopathology, Mayo Clinic, Rochester, MN

Peripheral blood flow cytometry is the test of choice to detect circulating neoplastic T-cells (Sézary cells) in patients with mycosis fungoides (MF) or Sézary syndrome (SS).<sup>1-4</sup> These 2 clinically distinct subtypes of cutaneous T-cell lymphoma

(CTCL) harbor neoplastic CD4<sup>+</sup> T cells with largely identical immunophenotypic properties, and share a staging scheme that includes blood ratings based on absolute numbers of Sézary cells (B0: <250, B1: ≥250 and <1000, and B2: ≥1000

cells/ $\mu\text{L}$ ).<sup>5</sup> Flow cytometric identification of Sézary cells also plays a critical role in distinguishing SS from reactive erythroderma, and in the assessment of therapy response. Unfortunately, the laboratory identification of Sézary cells has remained challenging, given that their most frequent immunophenotypic characteristics (lack of CD7 and CD26 expression) are shared with common benign CD4<sup>+</sup> T-cell subsets, whereas immunophenotypic properties specific for neoplasia (eg, diminished CD3 expression) are typically subtle or not identified.


We recently reported a novel flow cytometry strategy to rapidly confirm clonality on gated CD4<sup>+</sup> T-cell subsets suspicious for Sézary cells<sup>6,7</sup> using a single antibody against 1 of 2 mutually exclusive T-cell receptor constant  $\beta$  chains (TRBCs).<sup>8</sup> We hereby present a semiautomated analysis to detect Sézary cells using an unsupervised clustering tool from the EuroFlow Consortium-sponsored analysis software (Infinicyt version 2.0; Cytognos, Salamanca, Spain)<sup>9</sup> to automatically segment all CD4<sup>+</sup> T cells into immunophenotypically distinct clusters and test each cluster for clonality based on TRBC1 expression. In a large cohort of clinical specimens, we show the superior test performance of this semiautomated approach compared to traditional gating strategies, and objectively report the complex immunophenotypic properties of Sézary cells as identified without a preconceived bias of an expected immunophenotype.

Peripheral blood from 136 patients with inflammatory dermatoses (143 samples), 191 patients with CTCL (343 samples), and 22 healthy donors (22 samples) was evaluated by flow cytometry at Mayo Clinic, Rochester, MN, using a single-tube Sézary cell panel (CD2, CD3, CD4, CD5, CD7, CD8, CD26, CD45, and TRBC1 [clone JOVI.1]) on a FACSCanto II flow cytometer (BD Biosciences), as previously described.<sup>6</sup> Forty-five patients with inflammatory dermatoses had reactive erythroderma (33%); and 178 specimens from patients with CTCL were collected after systemic therapy or extracorporeal photopheresis (52%). Patients with MF were studied at stage I (125 samples), stage II (41 samples), stage III (30 samples), and stage IV (35 samples) of disease. Semiautomated analysis was performed in Infinicyt, including manual gating of CD4<sup>+</sup>/CD8<sup>-</sup> T cells per previously described strategies,<sup>3</sup> automated density-based clustering of CD4<sup>+</sup> T cells, automated detection of Sézary cell clusters exhibiting <15% or >85% TRBC1-positive cells, visual identification of Sézary cell clusters with dim TRBC1 expression, and expert gating upon suspicion of missed Sézary cells (Supplemental Figure 1, available on the *Blood* website). Absolute counts were calculated based on the absolute lymphocyte count obtained on a Sysmex XN-350 analyzer (Kobe). This study was approved by the Mayo Clinic Institutional Review Board.

Sézary cells were identified in 170 of 343 (50%) samples from patients with CTCL (range, 11-32 851 cells/ $\mu\text{L}$ ). Small TRBC-restricted CD4<sup>+</sup> T-cell subsets were also detected in a minority of benign samples (7 of 143 patients with inflammatory dermatoses; 4.9%; range, 11-155 cells/ $\mu\text{L}$ ), a finding not associated with erythroderma ( $P = .28$ ,  $\chi^2$ ) and consistent with T-cell clones of uncertain significance<sup>6,10</sup> (Figure 1A-D). Of 177 positive samples, 88% were automatically detected based on <15% or >85% TRBC1 positivity, 4% were visually identified based on dim TRBC1 expression, and 8% required expert gating of very small Sézary cell populations (range, 11-155 cells/ $\mu\text{L}$ ). When

compared to the gating strategy proposed by the European Organization for Research and Treatment of Cancer<sup>11</sup> (EORTC; maximum of CD4<sup>+</sup>/CD7<sup>-</sup> and CD4<sup>+</sup>/CD26<sup>-</sup> T cells), our approach produced very similar Sézary counts in the B2 high tumor burden range but resulted in significant discrepancies in the B0-B1 range (Figure 1E). Notably, 18 of 62 samples (29%) rated B1 by EORTC gating lacked detectable clonal T cells by TRBC1 analysis, and 50 of 205 (24%) samples rated B0 by EORTC gating had small CD4<sup>+</sup> T-cell clones by TRBC1 (Figure 1F). In samples from patients with inflammatory dermatoses, our semiautomated approach produced lower false Sézary cell counts (median, 0 vs 92 cells/ $\mu\text{L}$ ;  $P < .0001$ , Mann Whitney test), more accurately classified samples as B0 (diagnostic specificity, 100% vs 92%;  $P = .0015$ , McNemar test), and confidently identified lower Sézary cell counts (limit of blank, 9 vs 345 cells/ $\mu\text{L}$ ) compared to EORTC gating (Figure 1G). Moreover, EORTC Sézary cell counts on CTCL samples without detectable TRBC-restricted T cells closely mirrored those of patients with inflammatory dermatoses (Supplemental Figure 2A) consistent with poor Sézary cell discrimination by the EORTC method. Thus, our semiautomated approach to identify TRBC-restricted Sézary cells successfully overcame the limited specificity of the EORTC method and improved the confident identification of Sézary cells in the B0-B1 range.


We then characterized the immunophenotypic properties of TRBC-restricted Sézary cells. Combined density maps of TRBC-restricted Sézary cells (Infinicyt) showed extensive overlap with benign CD4<sup>+</sup> T cells from patients with inflammatory dermatoses on 2-dimensional projections of routinely evaluated antigens (Figure 2A), and on commonly used dimensionality reduction plots (Supplemental Figure 2B-C). Remarkably, there were no salient discriminatory features on the commonly used CD7 vs CD26 projection. The most distinct property of Sézary cells was variable loss of CD3 and CD2 expression, which was often but not always outside 2 standard deviations of the density contours of benign CD4<sup>+</sup> T cells. Unsupervised k-means metaclustering (Qlucore Omics Explorer version 3.7, Qlucore) revealed up to 5 immunophenotypically distinct classes of Sézary cells containing at least 5 clusters each, with an additional class of interest (class 6) identified on supervised analysis: class 1 (CD7<sup>-</sup>/CD26<sup>-</sup>), class 2 (CD7<sup>+</sup>/CD26<sup>-</sup>), class 3 (CD2<sup>dim/-</sup>), class 4 (CD3<sup>dim/-</sup>), class 5 (CD7<sup>-</sup>/CD26<sup>+</sup>), and class 6 (CD7<sup>+</sup>/CD26<sup>+</sup>) (Figure 2B). Slightly dim CD3 expression was a feature commonly encountered in all classes, except for class 4 where CD3 was overtly decreased to absent. Although most (59%) samples had only 1 class of Sézary cells, the presence of multiple (2, 3, or 4) classes in a single specimen was a frequent finding (29%, 8%, and 4% of positive specimens, respectively) (Figure 2C). Classes 2 and 4, when present, tended to be more abundant than the prototypic class 1 Sézary cell (adjusted  $P < .05$  for both); and classes 5 and 6 usually comprised a smaller fraction of total Sézary cells (adjusted  $P < .0001$  and  $P < .01$ , respectively) (Kruskal-Wallis and Dunn tests) (Figure 2D). Interestingly, class 2 was more frequently encountered in SS compared to MF ( $P < .05$ ), and class 4 was exclusively observed in MF (not statistically significant) (Fisher exact tests) (Figure 2E). Of 29 CTCL patients with more than 1 positive sample (average, 4.4 samples per patient; range, 2-11 samples), 15 (52%) showed a change in Sézary cell classes comprising the clonal subset, whereas shifts in the relative and absolute



**Figure 1. Semiautomated detection of constant beta chain-restricted CD4<sup>+</sup> T-cell clusters rapidly identifies Sézary cells, with improved specificity over traditional gating strategies.** (A) Blood samples from healthy donors (left), patients with benign inflammatory dermatoses (middle), and patients with CTCL (right) were analyzed by flow cytometry using automated clustering and calculation of percentage TRBC1-positive events per cluster. Individual CD4<sup>+</sup> T-cell clusters (dots) were classified as clonal (red) or nonclonal (blue) based on pre-established clonality thresholds (dotted lines) or monophasic TRBC1-dim expression on visual analysis. (B) Absolute numbers of Sézary cells (clonal CD4<sup>+</sup> T cells) per microliter in patients with CTCL. (C) Results of blood ratings based on total clonal CD4<sup>+</sup> T cells. (D) Number of clonal CD4<sup>+</sup> T-cell clusters identified per sample on patients with CTCL. (E) Absolute numbers of Sézary cells estimated by clustering and TRBC1 expression compared to the gating strategy recommended by the EORTC. (F) Blood ratings estimated by the EORTC-recommended gating strategy compared to automated clustering and TRBC1 expression. (G) Diagnostic specificity for B1-B2 blood rating (%B0) and analytical specificity (limit of blank) of the EORTC-recommended gating strategy (left) compared to automated clustering and TRBC1 staining (right); based on patients with inflammatory dermatoses.

abundance of each Sézary cell class were frequently encountered (Figure 2F). Overall, these findings highlight the complex and frequently unstable immunophenotypic properties of Sézary cells which often overlap with those of reactive CD4<sup>+</sup> T cells, making them challenging to identify by comprehensive immunophenotyping alone without clonality assessment.

In conclusion, we show the optimal test performance of a semiautomated flow cytometry strategy for objectively identifying Sézary cells based on the detection of immunophenotypically distinct CD4<sup>+</sup> T-cell clusters exhibiting a clonal TRBC1 expression pattern. This approach results in improved test specificity and confident Sézary cell count reporting in the



**Figure 2.**

**Figure 2. Immunophenotypic features of constant beta chain–restricted Sézary cells showing complex multicenter composition and overlap with benign T-cell subsets.** (A) Density plots showing extensive immunophenotypic overlap between constant beta chain–restricted Sézary cells from patients with CTCL (red lines), and benign CD4<sup>+</sup> T cells from patients with inflammatory dermatoses (blue lines). Double lines depict 2 standard deviations. (B) Metaclustering (K-means) of constant beta chain–restricted Sézary cell clusters showing distinct immunophenotypic subtypes on 2-dimensional projections (described as class 1-6). (C) Circos plot depicting the incidence (segments) and observed combinations (ribbons) of Sézary cell metaclusters. (D) Contribution of each immunophenotypic class of Sézary cell to the total Sézary cell count based on absolute numbers of cells per microliter (left) and percentage of total Sézary cells (right). Box and whiskers: median, interquartile range and total range. Asterisks (\*), (\*\*), and (\*\*\*) indicate  $P < .05$ ,  $P < .01$ , and  $P < .0001$ , respectively, compared to class 1. (E) Incidence of each immunophenotypic class of Sézary cell (first positive sample only) in patients with mycosis fungoides (grey bars) and Sézary syndrome (black bars). Asterisk (\*) indicates  $P < .05$ . (F) Examples of patients with Sézary syndrome (top) and mycosis fungoides (bottom), showing variability of Sézary cell metacluster composition and abundance on serial samples.

B0-B1 range over the proposed EORTC strategy. In addition, we describe the immunophenotypic spectrum of Sézary cells based on the unbiased identification of clonal CD4<sup>+</sup> T cells in patients with CTCL, including demonstration of the extensive immunophenotypic overlap with reactive T cells and a complex mixture of immunophenotypically distinct neoplastic subsets. Limitations of this strategy include inability to detect very rare circulating CD8<sup>+</sup> and CD4<sup>+</sup>/CD8<sup>+</sup> CTCLs (excluded per gating design), and subjectivity in discerning small CD4<sup>+</sup> T-cell clones of uncertain significance from low-level/minimal involvement by CTCL. Overall, our results support the upfront and standardized assessment of TRBC1 expression in the routine flow cytometric identification of Sézary cells.

## Authorship

Contribution: P.H., H.O., and G.O. designed the study; P.H. and G.O. performed case search and flow cytometric analysis; P.H. analyzed the data and wrote the manuscript; J.S., G.O., and H.O. contributed to the data analysis design; and D.J., M.S., and J.S. edited the manuscript and provided project support and expert knowledge.

Conflict-of-interest disclosure: The authors declare no competing financial interests.

ORCID profile: H.O., [0000-0001-5783-4485](https://orcid.org/0000-0001-5783-4485).

Correspondence: Pedro Horna, Mayo Clinic, 200 First Street SW, Rochester, MN 55905; email: [horna.pedro@mayo.edu](mailto:horna.pedro@mayo.edu).

## Footnotes

Submitted 24 June 2022; accepted 21 August 2022; prepublished online on *Blood* First Edition 29 August 2022.

The online version of this article contains a data supplement.

## REFERENCES

1. Pulitzer MP, Horna P, Almeida J. Sézary syndrome and mycosis fungoides: an overview, including the role of immunophenotyping. *Cytometry B Clin Cytom*. 2021;100(2):132-138.

2. Najidh S, Tensen CP, van der Sluijs-Gelling AJ, et al. Improved Sézary cell detection and novel insights into immunophenotypic and molecular heterogeneity in Sézary syndrome. *Blood*. 2021;138(24):2539-2554.
3. Horna P, Wang SA, Wolniak KL, et al. Flow cytometric evaluation of peripheral blood for suspected Sézary syndrome or mycosis fungoides: international guidelines for assay characteristics. *Cytometry B Clin Cytom*. 2021;100(2):142-155.
4. Vermeer MH, Moins-Teisserenc H, Bagot M, Quaglino P, Whittaker S. Flow cytometry for the assessment of blood tumour burden in cutaneous T-cell lymphoma: towards a standardized approach. *Br J Dermatol*. 2022;187(1):21-28.
5. Hristov AC, Tejasvi T, W RA. Cutaneous T-cell lymphomas: 2021 update on diagnosis, risk-stratification, and management. *Am J Hematol*. 2021;96(10):1313-1328.
6. Horna P, Shi M, Jevremovic D, Craig FE, Comfere NI, Olteanu H. Utility of TRBC1 expression in the diagnosis of peripheral blood involvement by cutaneous T-cell lymphoma. *J Invest Dermatol*. 2021;141(4):821-829 e822.
7. Shi M, Jevremovic D, Otteson GE, Timm MM, Olteanu H, Horna P. Single antibody detection of T-cell receptor alphabeta clonality by flow cytometry rapidly identifies mature T-cell neoplasms and monotypic small CD8-positive subsets of uncertain significance. *Cytometry B Clin Cytom*. 2020;98(1):99-107.
8. Maciocia PM, Wawrzyniecka PA, Philip B, et al. Targeting the T cell receptor beta-chain constant region for immunotherapy of T cell malignancies. *Nat Med*. 2017;23(12):1416-1423.
9. Pedreira CE, Costa ES, Lecrevisse Q, van Dongen JJM, Orfao A. Overview of clinical flow cytometry data analysis: recent advances and future challenges. *Trends Biotechnol*. 2013;31(7):415-425.
10. Shi M, Olteanu H, Jevremovic D, et al. T-cell clones of uncertain significance are highly prevalent and show close resemblance to T-cell large granular lymphocytic leukemia. Implications for laboratory diagnostics. *Mod Pathol*. 2020;33(10):2046-2057.
11. Scarisbrick JJ, Hodak E, Bagot M, et al. Blood classification and blood response criteria in mycosis fungoides and Sézary syndrome using flow cytometry: recommendations from the EORTC cutaneous lymphoma task force. *Eur J Cancer*. 2018;93:47-56.

<https://doi.org/10.1182/blood.2022017548>

© 2022 by The American Society of Hematology