

Check for updates

Direct oral anticoagulant (DOAC) interference in hemostasis assays

Karen A. Moser¹ and Kristi J. Smock²

¹Department of Pathology, University of Utah School of Medicine, Salt Lake City, UT; and ²ARUP Institute for Clinical and Experimental Pathology, Salt Lake City, UT

Direct oral anticoagulants (DOACs) are a group of direct coagulation factor inhibitors including both direct thrombin inhibitors and direct factor Xa inhibitors. These medications may cause hemostasis assay interference by falsely increasing or decreasing measured values, depending on the analyte. Considering the potential for DOAC interference in a variety of hemostasis assays is essential to avoid erroneous interpretation of results. Preanalytic strategies to avoid DOAC interference include selecting alternatives to clot-based hemostasis assays in patients taking DOACs when possible and sample collection timed when the patient is off anticoagulant therapy or at the expected drug trough. Clinical laboratories may also provide educational materials that clearly describe possible interferences from DOAC, develop testing algorithms to aid in detection of DOAC in submitted samples, use DOAC-neutralizing agents to remove DOACs before continuing with testing, and write interpretive comments that explain the effects of DOAC interference in hemostasis tests. Using a combination of the described strategies will aid physicians and laboratorians in correctly interpreting hemostasis and thrombosis laboratory tests in the presence of DOACs.

LEARNING OBJECTIVES

- Describe the patterns of interference in hemostasis assays due to direct thrombin inhibitors and direct Xa inhibitors
- List potential strategies physicians and clinical laboratories can use to decrease DOAC interference in hemostasis
 assays

Since the introduction of the direct oral anticoagulants (DOACs) in the 2010s, clinical laboratories have struggled to mitigate the effects of these drugs in clot-based hemostasis and thrombosis assays. DOACs are a group of direct coagulation factor inhibitors that include both direct thrombin inhibitors (dabigatran) and direct Xa inhibitors (rivaroxaban, apixaban, edoxaban).¹ The DOACs may cause assay interference by falsely increasing or decreasing measured values, depending on the analyte. Data from both individual laboratory studies and external quality assessment programs describe expected patterns with a variety of reagents (Table 1).¹⁻⁸ embolism, prompting testing for lupus anticoagulant (LA) and treatment with rivaroxaban. Laboratory results for both her initial LA profile as well as follow-up testing 12 weeks later are included in Table 2.

Case 2: A 37-year-old man sought treatment for an unprovoked pulmonary embolism and was evaluated for thrombophilia risk factors, including LA, while taking apixaban therapy. Laboratory results for his initial LA panel are included in Table 2. An anti-Xa activity assay calibrated for unfractionated heparin (UFH) showed measurable anti-Xa activity in the plasma sample submitted for the LA profile.

Patterns of DOAC interference in hemostasis assays

For basic hemostasis assays such as prothrombin time (PT) and activated partial thromboplastin time (aPTT), response to DOACs varies considerably by drug, drug concentration, and reagent. In general, the aPTT may be prolonged with dabigatran but does not tend to be prolonged by direct Xa inhibitors.^{12,9} A normal aPTT is insufficient to exclude the

CLINICAL CASES

Case 1: A 54-year-old woman was diagnosed with pneumonia and was bedridden due to the severity of her illness. A few weeks later, she developed a pulmonary

Table 1	. Patterns	of DOAC	interference	in hemostasis,	thrombosis /	assays
---------	------------	---------	--------------	----------------	--------------	--------

Expected change	Assays	Notes
Clotting time prolongation	• aPTT (dabigatran > direct Xa inhibitors) • PT (rivaroxaban > edoxaban > apixaban) • Thrombin time (dabigatran)	Effects on clotting times are reagent dependent. aPTT and PT mixing tests are expected to show incom- plete correction in the presence of DOACs.
False increase	 Clot-based protein C activity Clot-based protein S activity Antithrombin activity (in factor IIa-based assays with dabigatran, in factor Xa-based assays with direct Xa inhibitors) Activated protein C resistance ratio 	 False increase in protein C, protein S, and antithrombin activities may result in misdiagnosis of a patient with true deficiency as normal. Falsely elevated activated protein C resistance ratio may result in misdiagnosis of a patient with factor V Leiden mutation as normal.
False decrease	 aPTT-based factor assays (VIII, IX, XI, XII) PT-based factor assays (II, V, VII, X) 	Dilutions in factor assays may show nonspecific inhibitor effect.
False positive (or potentially false negative)	• LA assays	Includes aPTT- and DRVVT-based assays, among other clotting time-based LA assays; effects are drug and reagent dependent.
No change	 Clauss fibrinogen activity (for most reagents, rare methods show false decrease in presence of high concentrations of dabigatran) D-dimer Chromogenic protein C activity Free and total protein S antigen Anticardiolipin, anti-β2GP1 ELISAs von Willebrand activity and antigen assays DNA-based assays (eg, factor V Leiden mutation, prothrombin G20210A mutation) 	

presence of dabigatran.¹ The PT may be prolonged with rivaroxaban and edoxaban but tends to be relatively insensitive to the presence of apixaban in all reagents evaluated to date.^{2,9} The difference between the observed PT prolongation in the clinical cases highlights the differential effects of rivaroxaban and apixaban on the PT. Dabigatran does not tend to prolong the PT.^{2,9} In mixing tests for both aPTT and PT, DOACs will appear as nonspecific inhibitors, with the inhibitor effect most pronounced at higher drug concentrations.¹ Fibrinogen activity tends to be unaffected by DOACs, particularly the commonly used Clauss fibrinogen assay, which uses plasma dilution as well as a high concentration of thrombin to convert patient fibrinogen to fibrin.2-6 Factitiously decreased fibrinogen activity was reported by some laboratories at higher dabigatran concentrations (385 and 744ng/mL) in an external quality assessment program.³ D-dimer assays are typically immunoassays, such as enzyme-linked immunosorbent assay (ELISAs) or latex immunoassays; these methods are not affected by the presence of DOACs.¹

Among specialized hemostasis assays, those evaluating for thrombophilia risk factors are at particular risk for interference, given that many of these assays are clot based and may be measured in patients who have thrombosis and are treated with anticoagulant therapy.¹⁰⁻¹² The clinical cases illustrate two typical examples of how LA assays are affected by DOACs. Dilute Russell viper venom time (DRVVT) and aPTT-based LA assays (eg, platelet neutralization procedure, hexagonal phospholipid neutralization) may show false-positive results in the presence of DOACs.^{10,12,13} False-positive DRVVT-based testing appears to be a particular risk with rivaroxaban, but there also appears to be risk of DRVVT false negatives for LA in samples containing apixaban.¹⁴ Just as with screening aPTT and PT, the mixing test steps in DRVVT and aPTT-based LA assays may show a nonspecific inhibitor pattern.¹³ Assays for clot-based protein C, protein

S, and antithrombin may be falsely increased in the presence of DOACs, which may result in a false-negative result.^{2-6,11} That is, a patient with a true deficiency may have a normal result if plasma is tested for protein C, protein S, or antithrombin activity using clot-based assays in the presence of DOACs. DOACs do not affect chromogenic protein C activity assays or total protein C antigen assays.^{1,11} Protein S activity assays are clot based and experience interference by DOACs, whereas free and total protein S antigen assays are immunoassays and do not show DOAC interference.^{3,4,15} Antithrombin activity assays may be either factor IIa based or factor Xa based; the design determines which DOACs will interfere. Factor IIa-based assays show false elevation with direct thrombin inhibitors, whereas factor Xa-based assays show false elevation with direct Xa inhibitors.³⁻⁶ Activated protein C resistance assays also have the potential to give falsenegative results in the presence of DOACs.^{11,16}

Other clot-based hemostasis assays, such as factor assays, may show decreased activity in the presence of DOACs.^{1,3,4,11} Hemostasis assays based on ELISA or other immunoassay methods (eg, von Willebrand factor antigen, solid-phase antiphospholipid antibodies) will not show interference by DOACs; the results of these assays will not be affected if performed in a sample from a patient receiving DOAC therapy.¹ Likewise, DNA-based assays (eg, factor V Leiden mutation, prothrombin G20210A mutation) will be unaffected by the presence of DOACs.

One particularly challenging area where direct Xa inhibitor interference can limit therapeutic monitoring of another drug is the special case of patients switching from a direct Xa inhibitor to therapy with UFH or low molecular weight heparin (LMWH).¹⁷ Early reports of efficacy of DOAC-Stop (D-S; Haematex Research) for removing the measurable anti-Xa activity effect of rivaroxaban and apixaban but not heparins raise the possibility that plasmas in patients transitioning between a direct Xa inhibitor and UFH

Test	Case 1: initial LA panel	Case 1: LA panel 12 weeks later	Case 2: initial LA panel	Reference interval
PT(s)	23.0	12.8	14.3	12.0-15.5
DRVVT screen(s)	57	35	91	33-44
DRVVT 1:1 mix(s)	52	NA	80	33-44
DRVVT confirm(s)	Negative	NA	Positive	Negative
aPTT screen(s)	61	38	119	32-48
TT(s)	15.7	NA	15.1	14.7-19.5
aPTT 1:1 mix(s)	49	NA	85	32-48
PNP	Negative	NA	Positive	Negative
Hex phos neutralization	Positive	NA	NA	Negative

DRVVT confirm, PNP, and hex phos neutralization are LA confirmatory reagents containing high phospholipid concentrations.

hex phos neutralization, hexagonal phase phospholipid neutralization test; NA, not applicable; PNP, platelet neutralization procedure; TT, thrombin time.

or LMWH could be treated with D-S prior to measurement with a UFH- or LMWH-calibrated anti-Xa activity assay or aPTT for therapeutic monitoring.¹⁷⁻¹⁹ It is important to note that the degree to which activated carbon adsorbs LMWH is incompletely understood, and further study in this area is needed before routinely using DOAC-neutralizing compounds in samples also containing LMWH.¹⁹

What can physicians do to avoid DOAC interference in hemostasis assays?

Physicians can choose from a few different practical strategies to minimize DOAC interference in hemostasis assays. The simplest and best way to avoid DOAC interference is to avoid ordering and performing hemostasis assays in patients taking DOACs. Another option would be to stop DOAC therapy for 2 to 3 days prior to collecting a sample for hemostasis testing; however, given the risk of thrombosis with interruption or change of treatment, this option may not be clinically feasible in most cases.^{11,12,20} Briefly transitioning to an anticoagulant with less assay interference, such as LMWH, could also be an option but also may not be practical or feasible in many cases. In some clinical situations, testing in patients taking DOACs may be desired (eg, in patients requiring indefinite anticoagulation or detecting an LA in the setting of unprovoked thrombosis to gain information about recurrence risk).¹² If testing is undertaken, physicians can consider both choice of assay and timing of sample collection to decrease the possibility of DOAC interference. For analytes in which there is a choice between a functional assay and an assay such as a chromogenic, an antigenic, or a DNA-based assay where interference is not expected, physicians can choose the assay without expected interference. For example, using this principle would favor mutation analysis for factor V Leiden over functional testing with activated protein C resistance assays in a patient taking DOACs in whom thrombophilia testing is desired. Simply substituting another test does not work in all cases in which alternate assays are available. Some tests that lack interference (such as protein C or

antithrombin antigen) do not detect rare deficiencies due to protein dysfunction and would not represent appropriate substitutes for functional assays. If a measurement must be made, testing at expected drug trough concentration is recommended: either 12 hours postdose for DOACs dosed twice daily or 24 hours postdose for DOACs administered once daily.^{9,14} LA test systems (aPTT and DRVVT based) have demonstrated interference from DOACs even at low drug concentrations; therefore, testing a sample collected at an expected drug trough concentration may not completely eliminate the possibility of interference.²¹

What can clinical laboratories do to avoid DOAC interference in hemostasis assays?

Clinical laboratories likewise have multiple options for mitigating the effects of DOAC interference in hemostasis and thrombosis assays.¹⁴ A laboratory can provide educational material to physicians and other health care providers using the laboratory's services that include clear descriptions of expected patterns of DOAC interference, similar to Table 1.22 Laboratories may also design testing algorithms that allow detection of interfering anticoagulants in plasma samples in cases where a patient's history of anticoagulant therapy is not available.^{11,13} In clinical case 1, the prolonged PT offered a clue to the presence of an anticoagulant medication, either warfarin or a direct Xa inhibitor. Detecting anti-Xa activity in the sample is another, more sensitive way to identify a direct Xa inhibitor in a plasma sample; thrombin time can be used to identify the presence of direct thrombin inhibitors and heparins.^{11,13} Laboratories should also provide interpretive result comments that explain the potential for DOAC interference in the assays performed, particularly for LA profiles as recommended in current International Society on Thrombosis and Haemostasis (ISTH) and Clinical and Laboratory Standards Institute (CLSI) guidelines.^{13,23} Some laboratories may elect not to complete testing for a LA in the event DOAC interference is detected and may simply issue a comment stating that DOAC interference renders results of LA testing uninterpretable.²³

Recently, adsorbing agents that can neutralize DOAC effects in plasma samples in the hemostasis laboratory have been developed, although these agents are not yet FDA approved. Clinical laboratories have extensive experience with removing heparin effects from plasma samples using heparinase or polybrene.24 Using DOAC-neutralizing agents is an attractive possibility that could be incorporated into laboratory workflows that already use heparin neutralizers in a similar fashion. Currently available DOAC-neutralizing agents in tablet form include DOAC-Stop (D-S; Haematex Research) and DOAC-Remove (D-R; 5-Diagnostics AG).¹¹ These agents are composed of activated carbon-containing adsorbent compounds. One D-S or D-R tablet must be added to 1mL patient plasma, with time allowed for adsorption of DOACs and subsequent centrifugation and removal of DOAC-free plasma for testing (Figure 1).^{10,18,25} DOAC-spiked plasmas that show falsepositive LA results generally convert to negative LA results with DRVVT- and aPTT-based testing following D-S treatment, with results comparable to neutralization with idarucizumab or andexanet alfa.10,19,26,27 One tablet of D-S in 1mL plasma can reportedly neutralize up to 708ng/mL apixaban, 1060ng/mL edoxaban, 1020ng/mL rivaroxaban, and 360ng/mL dabigatran.¹⁰ D-S or D-R has also been reported to decrease DOAC interference in chromogenic and clot-based factor VIII activity,18,26 1-stage factor IX activity,²⁶ thrombin generation assays,²⁸ activated protein C



Figure 1. DOAC-neutralizing processes with tablet-based adsorbing agents (A) and filter systems (B).

resistance,^{25,29} and antithrombin activity.^{29,30} D-S has been reported to cause decreased factor activities in treated plasmas, raising the possibility of adsorption of coagulation factors in addition to DOACs; however, D-S does not produce aPTT prolongation with all reagents tested.^{10,31} Caution in result interpretation is recommended when using DOAC-neutralizing agents in the clinical laboratory as complete neutralization was not observed in all cases.^{10,18} One even more recently described option, DOAC Filter (Diagnostica Stago), is a cartridge containing a solid phase designed to extract DOACs through noncovalent binding when 600µL citrated plasma is added and centrifuged at 300×g for 15 minutes.³² DOAC Filter treatment resulted in undetectable levels of DOACs as assessed by quantitative assays and no change in selected hemostasis assays (PT, aPTT, fibrinogen, antithrombin activity, protein C activity, activated protein C resistance, LA tests) in postfiltration samples.³² Further study is needed to better describe the efficacy of DOAC-neutralizing reagents and the risk for inadvertent adsorption of other coagulation factors from plasma in the clinical laboratory. It is also important to remember that performing additional assays to detect DOACs and use of neutralizing agents adds cost and time to testing and requires additional patient plasma. Laboratories will need to carefully consider how these methods would fit into existing test workflows.

CLINICAL CASE (Continued)

Case 1: This patient had an initial positive LA profile with a follow-up profile 12 weeks later in which an LA was not detected; both profiles were performed while the patient was receiving rivaroxaban. Differential diagnostic considerations for this pattern include a transient LA, differences in timing of draw (eg, first profile drawn at drug peak and second profile drawn at drug trough), or nonadherence to the rivaroxaban regimen at the time of the second profile. The prolonged PT in the first LA profile suggests one of the latter 2 possibilities.

Case 2: This case demonstrates the relative insensitivity of the PT to apixaban. The measurable direct Xa activity suggests that the positive LA results obtained were caused by apixaban interference. A repeat sample collected when the patient is not taking apixaban, repeat sample collected at expected drug trough, or treating the current sample with a DOAC-neutralizing agent and repeating the LA profile are possible options to decrease apixaban interference.

Conclusion

Considering the potential for DOAC interference in a variety of hemostasis assays is essential to avoid erroneous interpretation of results. Preanalytic strategies to avoid DOAC interference include avoiding clot-based hemostasis assays in patients taking DOACs and sample collection timed when the patient is not taking anticoagulant therapy or at the expected drug trough (as opposed to peak or random samples). Strategies clinical laboratories may employ to avoid DOAC interference include providing educational materials that clearly describe possible interferences from DOACs, developing testing algorithms to aid in the detection of DOACs in submitted samples, using DOAC-neutralizing agents to remove DOACs before continuing with testing, and writing interpretive comments that explain the effects of DOAC interference in hemostasis tests. Using a combination of the described strategies will aid physicians and laboratorians in correctly interpreting hemostasis and thrombosis laboratory tests in the presence of DOACs.

Conflict-of-interest disclosures

Karen A. Moser: disclose no relevant conflicts of interest. Kristi J. Smock: disclose no relevant conflicts of interest.

Off-label drug use

Karen A. Moser: nothing to disclose. Kristi J. Smock: nothing to disclose.

Correspondence

Karen A. Moser, ARUP Laboratories, 500 Chipeta Way, Mail Stop 115-G04, Salt Lake City, UT 84108; e-mail: karen.moser@hsc.utah .edu.

References

- Gosselin RC, Adcock DM, Bates SM, et al. International Council for Standardization in Haematology (ICSH) recommendations for laboratory measurement of direct oral anticoagulants. *Thromb Haemost*. 2018;118(3):437-450.
- Gosselin R, Grant RP, Adcock DM. Comparison of the effect of the anti-Xa direct oral anticoagulants apixaban, edoxaban, and rivaroxaban on coagulation assays. Int J Lab Hematol. 2016;38(5):505-513.
- Bonar R, Favaloro EJ, Mohammed S, Pasalic L, Sioufi J, Marsden K. The effect of dabigatran on haemostasis tests: a comprehensive assessment using in vitro and ex vivo samples. *Pathology*. 2015;47(4):355-364.
- Bonar R, Favaloro EJ, Mohammed S, et al. The effect of the direct factor Xa inhibitors apixaban and rivaroxaban on haemostasis tests: a comprehensive assessment using in vitro and ex vivo samples. *Pathology*. 2016;48(1):60-71.
- Van Blerk M, Bailleul E, Chatelain B, et al. Influence of apixaban on commonly used coagulation assays: results from the Belgian national external quality assessment scheme. *Int J Lab Hematol.* 2017;39(4):402-408.
- Van Blerk M, Bailleul E, Chatelain B, et al. Influence of dabigatran and rivaroxaban on routine coagulation assays: a nationwide Belgian survey. *Thromb Haemost.* 2015;113(1):154-164.
- Hillarp A, Gustafsson KM, Faxälv L, et al. Effects of the oral, direct factor Xa inhibitor apixaban on routine coagulation assays and anti-FXa assays. J Thromb Haemost. 2014;12(9):1545-1553.
- Hillarp A, Baghaei F, Fagerberg Blixter I, et al. Effects of the oral, direct factor Xa inhibitor rivaroxaban on commonly used coagulation assays. J Thromb Haemost. 2011;9(1):133-139.
- Douxfils J, Ageno W, Samama C-M, et al. Laboratory testing in patients treated with direct oral anticoagulants: a practical guide for clinicians. J Thromb Haemost. 2018;16(2):209-219.
- De Kesel PM, Devreese KMJ. Direct oral anticoagulant adsorption: impact on lupus anticoagulant testing—review of the literature and evaluation on spiked and patient samples. J Thromb Haemost. 2020;18(8):2003-2017.
- Siriez R, Dogné J-M, Gosselin R, Laloy J, Mullier F, Douxfils J. Comprehensive review of the impact of direct oral anticoagulants on thrombophilia diagnostic tests: practical recommendations for the laboratory. *Int J Lab Hematol.* 2021;43(1):7-20.
- Tripodi A, Cohen H, Devreese KMJ. Lupus anticoagulant detection in anticoagulated patients: guidance from the Scientific and Standardization Committee for lupus anticoagulant/antiphospholipid antibodies of the International Society on Thrombosis and Haemostasis. J Thromb Haemost. 2020;18(7):1569-1575.
- Devreese KMJ, de Groot PG, de Laat B, et al. Guidance from the Scientific and Standardization Committee for lupus anticoagulant/antiphospholipid antibodies of the International Society on Thrombosis and Haemostasis:

update of the guidelines for lupus anticoagulant detection and interpretation. J Thromb Haemost. 2020;18(11):2828-2839.

- Favaloro EJ, Mohammed S, Curnow J, Pasalic L. Laboratory testing for lupus anticoagulant (LA) in patients taking direct oral anticoagulants (DOACs): potential for false positives and false negatives. *Pathology*. 2019;51(3):292-300.
- Smock KJ, Plumhoff EA, Meijer P, et al. Protein S testing in patients with protein S deficiency, factor V Leiden, and rivaroxaban by North American Specialized Coagulation Laboratories. *Thromb Haemost*. 2016;116(1):50-57.
- Moore GW, Van Cott EM, Cutler JA, Mitchell MJ, Adcock DM; Subcommittee on Plasma Coagulation Inhibitors. Recommendations for clinical laboratory testing of activated protein C resistance; communication from the SSC of the ISTH. J Thromb Haemost. 2019;17(9):1555-1561.
- Douxfils J, Adcock DM, Bates SM, et al. 2021 Update of the International Council for Standardization in Haematology recommendations for laboratory measurement of direct oral anticoagulants. *Thromb Haemost*. 2021;121(8):1008-1020.
- Platton S, Hunt C. Influence of DOAC Stop on coagulation assays in samples from patients on rivaroxaban or apixaban. Int J Lab Hematol. 2019;41(2):227-233.
- Frans G, Meeus P, Bailleul E. Resolving DOAC interference on aPTT, PT, and lupus anticoagulant testing by the use of activated carbon. *J Thromb Hae*most. 2019;17(8):1354-1362.
- Godier A, Dincq A-S, Martin A-C, et al. Predictors of pre-procedural concentrations of direct oral anticoagulants: a prospective multicentre study. *Eur Heart J.* 2017;38(31):2431-2439.
- 21. Ratzinger F, Lang M, Belik S, et al. Lupus-anticoagulant testing at NOAC trough levels. *Thromb Haemost*. 2016;116(2):235-240.
- 22. College of American Pathologists. Hematology and Coagulation Checklist. College of American Pathologists; 2021.
- CLSI. Laboratory Testing for the Lupus Anticoagulant; Approved Guideline. CLSI document H60-A. Clinical and Laboratory Standards Institute; 2014.
- 24. Harenberg J, Reichel T, Malsch R, Hirsh J, Rustagi P. Multicentric evaluation of heparinase on aPTT, thrombin clotting time and a new PT reagent based on recombinant human tissue factor. *Blood Coagul Fibrinolysis*. 1996;7(4):453-458.
- Jilma-Stohlawetz P, Lysy K, Sunder-Plassmann R, et al. Limitations of a calibrated, quantitative APC-R assay under routine conditions. *Int J Lab Hema*tol. 2021;43(2):318-323.
- 26. Favaloro EJ, Gilmore G, Bonar R, et al. Reducing the effect of DOAC interference in laboratory testing for factor VIII and factor IX: a comparative study using DOAC Stop and andexanet alfa to neutralize rivaroxaban effects. *Hαemophilia*. 2020;26(2):354-362.
- Baker SA, Jin J, Pfaffroth C, Vu T, Zehnder JL. DOAC-Stop in lupus anticoagulant testing: direct oral anticoagulant interference removed in most samples. *Res Pr Thromb Haemost*. 2021;5(2):314-325.
- Kopatz WF, Brinkman HJM, Meijers JCM. Use of DOAC Stop for elimination of anticoagulants in the thrombin generation assay. *Thromb Res.* 2018;170:97-101.
- Favresse J, Lardinois B, Sabor L, et al. Evaluation of the DOAC-Stop* procedure to overcome the effect of DOACs on several thrombophilia screening tests. *TH Open*. 2018;2(2):e202-e209.
- Ząbczyk M, Natorska J, Kopytek M, Malinowski KP, Undas A. The effect of direct oral anticoagulants on antithrombin activity testing is abolished by DOAC-stop in venous thromboembolism patients. Arch Pathol Lab Med. 2021;145(1):99-104.
- Riva N, Vella K, Hickey K, et al. The effect of DOAC-Stop(*) on several oral and parenteral anticoagulants. Int J Lab Hematol. 2021;43(4):171-175.
- Sevenet PO, Cucini V, Hervé T, et al. Evaluation of DOAC Filter, a new device to remove direct oral anticoagulants from plasma samples. Int J Lab Hematol. 2020;42(5):636-642.

© 2021 by The American Society of Hematology DOI 10.1182/hematology.2021000241